
WARNINGS

1. Discontinue the use of the device if damage to the device is identified.
2. Use with care.
3. For single use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity and result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also contaminate the device and result in patient infection which may lead to patient illness or death. Manufacturer makes no performance claim for product that is reused, reprocessed, or resterilized.


SHAPE, STRUCTURE AND PRINCIPLES

1. Structural drawing

(1) Main body

(2) Accessory (Stopcock)

2. Summary of device

REF_GM-30N,GM-30NF,GM-30NS : INFLATION DEVICE with 310mm high-pressure extension tube length.
 REF_GM-30NL,GM-30NFL,GM-30NSL : INFLATION DEVICE with 766mm high-pressure extension tube length.
 Stopcock: GM30N, GM30NL = Three-way stopcock
 GM30NF, GM30NFL = Four-way stopcock

The device is a handheld inflation device equipped with a pressure gauge and is used to inflate and deflate balloon catheters by manipulating the plunger and the lever. The maximum allowable pressure is 30 atm/bar (441 psi). The accuracy of the pressure gauge is within $\pm 3\%$ of the maximum allowable infusion pressure. Stopcock is attached as an accessory.

INTENDED USE AND EFFICACY OR EFFECTS

This device is recommended for use with balloon dilation catheters to create and monitor pressure in the balloon, and to deflate the balloon.

CONTRAINDICATIONS AND PROHIBITIONS

- (1) Do not apply sudden positive or negative pressure.
- (2) Do not use with drugs containing organic solvents, fat emulsions and oily ingredients.

COMPLICATIONS

The following complications may occur as a consequence of use of this device:

- Vascular rupture, perforation, injury
- Death
- Infection
- Vascular dissection
- Embolism (air, tissue, thrombus)
- Bleeding complications
- Drug reaction, allergic reaction to contrast media, etc.

OPERATION OR USAGE METHOD

1. Usage method

NOTE: GM-30NL,GM-30NFL,GM-30NSL model, prior to use remove the clip that ties the high-pressure extension tube together.

- (1) The device should be carefully examined prior to use for defects to ensure proper performance.
- (2) Prepare a diluted contrast medium (mixed mainly with heparinized saline solution) that the balloon catheter manufacturer recommends.
- (3) Pull the plunger slowly while pressing the lever to aspirate an appropriate amount of the diluted contrast medium through the stopcock.
- (4) Hold the device upright so that the pressure gauge is above, turn the plunger clockwise slowly, and completely remove the air in the device.
- (5) Confirm that the amount of the diluted contrast medium recommended by the balloon catheter manufacturer (typically around 5 mL) is in the device.
- (6) Connect the stopcock with care so that the air may not get into the expansion port of the balloon catheter.
- (7) To inflate the balloon catheter, turn the plunger clockwise.
- (8) To deflate the balloon catheter, turn the plunger counterclockwise (the balloon catheter cannot be inflated with the lever depressed).

2. Precautions concerning usage method

- (1) Do not use the device in hospitals where urgent surgical intervention cannot be provided in case of complications that may lead to injury or severe life-threatening complications.
- (2) Use the device with the air in it completely removed. [Complications such as air embolism may occur.]
- (3) Regularly check for the damage of the device, loosening of the joint, and leakage of drug solution.
- (4) Copper ion may be eluted from the wetted part of the device, resulting in discoloration of the contrast medium into green, etc.

Storage method

Store the product in a room temperature location not exposed to high temperature and humidity or direct sunlight and take proper precautions to ensure product does not contact water.

NAME OF MANUFACTURER

Manufactured by : Goodman Co., Ltd.
 5F KDX Nagoya Sakae Building, 4-5-3
 Sakae, Naka-ku, Nagoya, Aichi 460-0008 Japan

Web : www.goodmankk.com/

Phone number : +81-52-269-5300

U.S.A. Local Contact : Infraredx, Inc.

28 Crosby Drive, Suite 100

Bedford, MA 01730

USA

Web : www.infraredx.com

Phone number : +1-800-596-3104

LABELING SYMBOL DEFINITIONS

	Catalogue number		Batch code
	Use-by date		Manufacturer
	Sterilized using ethylene oxide		
	Consult instructions for use		
	Do not re-use		Do not re-sterilize
	Do not use if package is damaged		
	Keep away from sunlight		Keep dry
	Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.		

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Fabricado por : Goodman Co., Ltd.