Goodtec Angiographic Diagnostic Catheter

ENGLISH

WARNING

- Patient selection should be considered given that use of this product may result in subacute thrombosis, cause intravascular complications or hemorrhagic complications.
- During procedure, perform administration of appropriate anticoagulant and antiplatelet therapy, as directed by lead physician.
- The product is a medical device and should only be used by physicians trained in the procedures of coronary angiography (CAG) and percutaneous transluminal coronary angioplasty (PTCA). The product should only be used in facilities capable of performing emergency coronary artery bypass grafting (CABG) as a provision against complications that may cause injury and serious complications that could prove to be life threatening.

Combination of medical devices

- Reference package insert for each medical device and pharmaceutical product used together with this product and use them with a full understanding of their characteristics to prevent damage to the product.
- Immediately cease use of any device or pharmaceutical product should any complications arise.

Method of use

- Select appropriate size and shape in consideration of the applicable case.
- All product insertion to target location and removal should be undertaken with an appropriately sized guide wire (this minimizes kinking and damage to product).
- Insertion of product through the sheath introducer should be undertaken with the guide wire inserted up to tip shape. For a pigtail product, the shape should be extended via guide wire after passing through the inducer (this eliminates damage to product shape or to catheter body).
- Care should be taken to ensure that the catheter tube tip does not cause damage to the vasculature during insertion.
- Intravascular manipulation of the product should be carefully performed under fluoroscopy, confirming both tip movement and positioning (kinking or damage may result from stress concentrated at the difference in hardness between the special soft middle tube and the adjoining tube for some types of models such as YUMIKO).
- Should any product resistance or a reduction in torque of the catheter be determined during operation, immediately cease usage and confirm cause (continued use may result in damage to vessel or the product).
- Immediately cease use and exchange product should kinking or damage to product be confirmed during use.
- Do not manipulate guide wire if kinking or damage to catheter is confirmed or guide wire is unable to advance (over manipulation may result in damage to catheter tip, body or vasculature).
- Before inserting any pharmaceutical products, confirm that there is no damage or blockage in product (damage may result due to excessive pressure).
- When removing product ensure that the guide wire has been inserted. If the guide wire is unable to be placed within product due to kinking or damage prior to removal, simultaneously remove product and sheath introducer.

CONTRAINDICATIONS FOR USE

- The product is sterilized, is not re-usable and must not be re-sterilized.

Inappropriate patients

- Patients with extensive peripheral disease which prevents the insertion of an appropriate size sheath (damage to vasculature may result).
- Patients in whom antiplatelet and/or anticoagulant therapy is contraindicated for reasons such as excessive prolongation of coagulation time (complications relating to hemorrhage may result).
- Patients with allergy to any pharmaceutical products used in the course of a PTCA procedure (complications relating to shock may develop).
- Patients who are pregnant or suspected of being pregnant (radiation exposure to fetus may occur).
- Patients with extremely low left ventricular function (outside guidelines for PTCA).
- Patients having gastrointestinal bleeding (further complications may arise).

Prohibition of use with medical devices

- Do not use any products comprising organic solvent such as ethanol, fat-based emulsions or oil-based contrast media (damage to the product or connector may result).

Prohibition of use

- Do not exceed the maximum pressure as indicated on the label when using an automated injection machine to inject contrast media or pharmaceutical product (damage to product may result).
- Do not apply heat, bend or make any additional holes to catheter tip (damage to product may result).

SHAPE / CONSTRUCTION


This product comprises of various catheter tubing attached to a shaped tip to form a disposable angiographic catheter. This product's tip is delivered to target location under fluoroscopic imaging to highlight the construction of the vasculature via a flushing of contrast media through the device. As a coronary sinus device, blood sampling is also performed.

1. Catheter tube section

The catheter tubing is comprised from a soft tip, middle tubing and shaft, and when necessary, soft tip, middle tubing and inducer. Further, the soft tip is radiopaque, and the shaft is made of a flexible, non-braided and braided section.


2. Proximal end

The proximal end comprises of a strain relief and connector.

Maximum Pressure Resistance: 7240kPa (1050 psi)
Catheter tube material: nylon

Major catheter types

PURPOSE, EFFICACY

The following three applications of the device are performed based on type of catheter selected:

1. Coronary catheterization is achieved by intravascular insertion of the catheter and injection of contrast media in order to undertake coronary angiography (CAG), left ventriculography, abdominal angiography, diastolic aortic pressure (OAP), left ventricular pressure (LVP), and cerebral angiography.
2. Coronary catheterization can be performed via one catheter for the purpose of obtaining coronary angiography, left ventriculography and left ventricular pressure.
3. Obtain blood specimens from different locations within the coronary artery and large vessel in order to undertake blood analysis to assist with diagnosis.

SPECIFICATIONS

1. Maintains tensile strength of the catheter tubing at 9.8N (1kgf) without any separation or cracking.
2. No connector separation at 9.8N (1kgf).
3. The connection between soft tip and mid-tubing maintains integrity at 2.9N (0.3kgf).
4. The recommended guide wire diameter for 1.40mm to 2.00mm (4Fr to 6Fr) is 0.97mm (0.038"), or depending on the type may be 0.89mm (0.035") and for 1.13mm (3.4Fr) is 0.64mm (0.025").

METHOD OF USE

The following is an explanation of general use. Further steps may be required or changed as per the experience of the treating physician.

1. Preparation prior to use

- Carefully remove the product from the packaging.
- Inspect and confirm the catheter tubing and hub connector to ensure no product deformity (exchange product should any abnormality be determined).
- Using a syringe containing heparinized saline flush the inner lumen of the catheter.

2. Insertion of product

(1) Angiographic use

- Securing the sheath introducer, insert an appropriately sized angiographic guide wire.
- Placing the guide wire through the product, insert catheter through the sheath and into the vasculature.
- Advance catheter under fluoroscopy to the target vessel.
- Remove guide wire once catheter has reached target location (always ensure guide wire is inserted in product when delivering pig-tail catheter) (this minimizes the incidence of kinking and damage).
- Insert contrast media or pharmaceutical product.
- When performing co-axial technique, remove the guide wire and carefully insert the micro catheter through the product (the micro catheter size should not be larger than the guide wire size which is compatible with the product).

(2) Coronary sinus use

- Inserting an appropriately sized guide wire through the product, insert into the vasculature via a sheath introducer.
- Insert product into the right atrium using the guide wire.
- Ensuring there is no slack in the shape of the product, torque in a right direction until tip is gently placed against the right atrial posterior wall.
- With the tip placed against the wall, slowly lower the body of the catheter, placing the tip against the entry of the coronary sinus.
- With the tip placed in the coronary sinus slowly lower the cath-

eter body, extending the shape and inserting the tip within the coronary sinus.

3. Removal of catheter

- Remove ensuring the guide wire is inserted through catheter to tip location.

PRECAUTIONS

Precautions prior to use

- Use product prior to date of expiration as noted on label.
- Confirm condition of product prior to use. If any damage or contamination is noted on product, immediately cease use and exchange product.
- When removing product from packaging, remove with the cardboard (damage to product or tip shape may otherwise occur).
- Immediately use product following removal from package.
- All use should be performed in a sterile environment.
- Ensure there is no slack between product and hub connection prior to use (leakage of contrast media may incur).
- Do not use excessive force when removing residual air through the connector (damage to connector may result).

Precautions during use

- Ensure product does not come into direct contact with scissors, surgical blades or needles during use. Further, should any leakage or abnormality of use be determined, immediately exchange the product.
- When using the inducer to place the device through a sheath introducer, secure the inducer by hand when inserting (not securing the inducer may result in inducer being inserted into the vasculature).
- Do not exceed maximum injection pressure during insertion of contrast media etc. Contrast media should be maintained at 37°C.
- Ensure that connector is appropriately connected.
- Careful product manipulation should be undertaken when there is potential for a lesion in the ostial of the vessel.

Defective product

Use of product may result in the following possible defective product occurrences:

- Kinking or breaking of shaft
- Detachment or burst of shaft
- Difficulty in removing catheter
- Difficulty in guide wire manipulation and operation
- Damage to connector

Adverse event

Possible adverse events include, but are not limited to, the following:

Complications associated with medicinal products such as contrast media, infection, hemorrhage/hematoma, artery embolism/occlusion, arterial dissection, arrhythmia including ventricular fibrillation, arterial damage/perforation, stroke, subintimal injection of contrast, arteriovenous fistula, pseudoaneurysm, hypotension (severe hypotension), blood loss from puncture site, coronary arterial dissection, hemorrhagic complications, acute myocardial infarction, aortic dissection, unstable angina, distal embolization, chills/fever, intravascular thrombosis, nausea/vomiting, behavioral disorder, kidney injury, obstruction, palpitation, tachycardia and bradycardia

Precautions after use

Dispose of this product as medical waste, and take measures to prevent possible spread of infection.

STORAGE METHOD, SHELF LIFE

1. Storage method

- Store the product in a room temperature location not exposed to high temperature and humidity or direct sunlight and take proper precautions to ensure product does not contact water.
- Avoid inclinations, vibrations and impacts (including during transportation) and store in a safe, stable environment.
- Do not store near chemicals or in areas where the device may be exposed to gases.

2. Shelf life

Use this product before the "expiration date" shown on the package label.

PACKAGING

1 to 5 units / box

PRODUCT WARRANTY DISCLAIMERS AND LIMITATION OF REMEDIES

WITH RESPECT TO ANY AND ALL GOODMAN PRODUCT(S) DEPICTED OR DESCRIBED IN THIS PUBLICATION, THERE ARE NO EXPRESS OR IMPLIED WARRANTIES OF ANY KIND WHATSOEVER, INCLUDING, AND WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. GOODMAN CO., LTD. AND ITS SUBSIDIARIES ("GOODMAN") SHALL NOT BE LIABLE FOR ANY DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES IN CONNECTION WITH THE USE, RE-USE, OR ANY OTHER ASPECT OF THE RELEVANT PRODUCT(S) EXCEPT AS EXPRESSLY PROVIDED UNDER APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND GOODMAN TO ANY REPRESENTATION OR WARRANTY. ANY PRODUCT OR SERVICE DESCRIPTION OR SPECIFICATIONS CONTAINED IN ANY GOODMAN PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE RELEVANT PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS OR OTHER WARRANTIES.

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Goodman Co., Ltd.