NIPROTM A2F

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INTENDED PURPOSE:

Preparation of dialysis fluid used in hemodialysis and related therapies.

INDICATIONS:

Treatment of patients with acute or chronic renal failure when hemodialysis is prescribed by a physician.

CLINICAL BENEFITS:

Sodium bicarbonate is a buffer for the treatment for acidosis. The use of powder sodium bicarbonate reduces the risk of microbial growth in the concentrate.

CONTRAINDICATIONS:

There are no special contraindications for the use of NiproCart for hemodialysis. Patients who are hypersensitive to sodium bicarbonate or any components used in the NiproCart should not be dialyzed with the product.

SIDE EFFECTS:

There are no side-effects reported for NiproCart.

STORAGE:

NiproCart must be stored at temperatures below 30°C, in a cool and dry place.

WARNINGS:

- To be used only by healthcare professionals. Home dialysis is excluded.
- Use only if the packaging is intact.
- Before use, please check the cartridge for powder leakage. Do not use the cartridge if you acknowledge any powder leakage through the bottom connector, or a significant powder leakage through the top connector after turning the cartridge upside down.
- Handle with care to prevent microbial contamination.
- Use immediately after opening.
- Bicarbonate concentrate solution should be used just after reconstitution of bicarbonate powder with dialysis water by the hemodialysis monitor.
- Conductivity of the bicarbonate concentrate can be guaranteed up to 24h in a dialysis session.
- Single-use product. Any unused product must be discarded. Using the product multiple times may lead to contamination of the product.
- Non-pyrogenic product.
- Do not use after the expiry date stated on the packaging.
- The deionized water used must be pyrogen-free and must meet EP and ISO23500-5 standards.
- For use with Nipro hemodialysis monitors and BBraun Dialog+, BBraun Dialog IQ, Baxter AK200, Baxter AK200S, Baxter AK98, Baxter Artis, Nikkiso DBB05, and Nikkiso BDD EXA.
- Connect the sodium bicarbonate cartridge with an appropriate acid concentrate, according to medical prescription.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is established.
- Dispose the product according to local and national regulations.

CHECK-UP AND MAINTENANCE:

Before a hemodialysis session, check that all connections (acid concentrate and cartridge) are correct. After each session, clean and disinfect the monitor according to the manufacturer's instructions.

MEDICAL INFORMATION

The use of the bicarbonate powder cartridge with the chosen acid concentrate (in the right proportions) will give the concentrations stated on the acid concentrate label. It is the physician's responsibility to choose the formula that each patient requires, to report any adverse effects that may arise during hemodialysis, and to ensure proper use of the product by the staff.

DILUTION AND CAPACITY

100 g of sodium bicarbonate powder will generate 1 L of sodium bicarbonate concentrate

ACID CONCENTRATE REQUIREMENTS:

Use with Nipro acid concentrates serial 300 and 500 or equivalent acid concentrate to be used with pure bicarbonate concentrate. Sodium bicarbonate concentrate must be diluted in the dilution ratio expressed in the acid hemodialysis bath indicated in the acid concentrate label.

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Nipro Australia Pty Ltd