## Guide Extension Catheter

### Instructions for Use

Please read these instructions carefully before use.

#### Indications

This device is inserted into coronary arteries and used to guide interventional devices, etc. to affected areas in coronary arteries during procedures such as percutaneous transluminal coronary angioplasty. Notably, this product is placed in front of affected areas that are difficult for interventional devices to reach, and then used to support an interventional device.

#### Contraindications

- Contraindicated patient populations
  - Patients with history of coronary spasm. [Acute coronary occlusion may occur.]
  - Use in a left main coronary trunk which is not being protected by a bypass or collateral circulation. [Acute coronary occlusion may occur.]
  - Patients contraindicated for antiplatelet therapy or anticoagulant therapy utilized when there is excessive prolonging of coagulation time. [It might become difficult to stop bleeding.]
  - Patients with severe allergies to any drugs needed for procedures, e.g. contrast agent. [An allergic reaction may be induced.]
  - Patients with an extremely low-functioning left ventricle. [Such patients have higher risks for death, myocardial infarction or ventricular fibrillation during and after surgery.]
- Instructions for use
  - Do not re-use, re-sterilization.
  - Do not use chemicals containing any organic solvents or oil-based contrast agents. [Deformation, deterioration or breakage may occur.]
  - When introducing the contrast media into this device, do not apply high-pressure infusion with an auto-injector or a similar tool. [High pressure may make this device protrude, resulting in damage to a vessel.]

#### Warning


- Instructions for use
  - This device must be used only at the institute where emergency coronary artery bypass surgery is conducted considering the possibility of occurrence of life-threatening and any other serious adverse events.

#### Shape, Structure

- Shape and Structure
 

This device is sterilized using ethylene oxide.

This device is comprised of a tip, distal shaft, proximal end protective shaft, and the proximal shaft which is made from metal wire.

- Radiopaque markers are on the proximal end of distal shaft.
- Depth markers are located 900 mm and 1000 mm from the tip.
- Compatible guide catheter inner diameter: 1.78 mm or more (0.070 inches or more)

#### Instructions for Use

- Points of confirmation before use
 

Before proceeding with treatment, carefully inspect all tools and equipment, including this device, and confirm they are suitably functional for the treatment to be performed.
- Preparation
  - Open the package, and take out the device together with the catheter cover.
  - Take out the device from the catheter cover.
  - Flush the inside of the distal shaft of this device using heparinized physiological saline. Also, moisten the surface of this device with heparinized physiological saline and confirm that the hydrophilic coating of the distal shaft is lubricated.
- Catheter Insertion
  - Insert the sheath introducer and guiding catheter by means of percutaneous puncture.
  - Insert the guidewire through the vascular catheter Y-connector connected to the guiding catheter, and advance it to the target area of the coronary artery.
  - Advance this device following the guidewire until it is in front of the hemostasis valve of the vascular catheter Y-connector connected to the guiding catheter.
  - Open the hemostasis valve and advance this device into the guiding catheter. When doing so, carefully advance this device once the depth marker reaches the location of the vascular catheter Y-connector.
  - Under fluoroscopy, advance this device beyond the distal tip of the guiding catheter until it reaches the target area within the coronary artery.
  - While holding this device, carefully advance the combined-use device (interventional device, etc.) following the guidewire until it reaches the affected part of the coronary artery. Then, perform the treatment for the affected area.
- Catheter Removal
  - Loosen the hemostasis valve of the vascular catheter Y-connector.
  - With the guidewire indwelling in the coronary artery, remove this device following the guidewire.

#### <Precautions on Instructions for use>

- Confirm before use that the size and configuration of this device are suitable for the conditions of the treatment. Confirm that the size and configuration of the combined-use device (device used in combination with this device) are suitable as well. [Blood vessels may be damaged, or this device or a combined-use device may be damaged.]
- In cases where there is only a limited gap between the combined-use device and the lumen of this device's distal shaft, to reduce the risk of an air embolism, open the hemostasis valve, confirm there is back-bleeding and then perform device delivery or pull-out.
- When the catheter is inside the body, always perform operations under fluoroscopy.
- Perform work inside blood vessels carefully, and stop operation even if a small amount of resistance is felt. Confirm the cause of the resistance under high-resolution fluoroscopy. Never proceed with an operation or withdraw the device until the cause of resistance is determined. [Blood vessels may be damaged, or device or guidewire breakage or damage may occur.]

- When using this device, take sufficient precautions to avoid damaging the walls of blood vessels.
- When inserting this device into the target blood vessel, be careful not to completely block the flow of blood.
- When inserting this device, do not apply any torque. [This may cause kinks, twists, or damage to the device.]
- For this device, the recommended depth of insertion into a coronary artery is more than 2 cm. [Insertion depths of 2 cm or less may not provide sufficient backup force.]
- Do not insert more than one combined-use device at a time. [This may break or damage this device or combined-use devices.]
- Do not use 2 or more of these devices at once in 1 guiding catheter.
- If intravascular pressure drops after this device is inserted, remove this device immediately.
- Do not advance this device 15 cm beyond the distal tip of the guiding catheter. [The device may catch onto the end of the guiding catheter or may become impossible to remove.]
- When inserting this device into an area where a stent has already been placed, proceed with insertion carefully so as not to damage or deform the stent. If this device gets caught on the lumen of the stent, do not push in or pull out this device with excessive force. [This may damage blood vessels, damage or displace the stent, or break or damage this device.]
- Do not withdraw an undeployed stent that has exited the tip of this device back into this device. To withdraw an undeployed stent, pull out both the stent and this device together back into the guiding catheter, and then remove. [The stent may catch onto the tip of this device, which can dislodge the stent.]
- When resistance is felt against a combined-use device inside this device or it gets stuck, do not pull out or push in the combined-use device or this device independent of one another. Always remove the combined-use device together with this device at the same time. [This may break or damage this device or combined-use devices.]

#### Caution

- Important basic cautions
  - Do not use this device in blood vessels with a diameter of less than 2.5 mm. [There is a risk of blood vessel damage, ischemia, or occlusion.]
  - Provide appropriate anticoagulant therapy while using this device.
  - If kinking or twisting occurs while operating this device, remove the whole system from the patient.
  - Do not use if package is damaged.
- Defects or adverse events
 

Use of this device carries the risk of the following defects or adverse events occurring.

  - Serious defects
    - Catheter damage
    - Catheter removal difficulty
  - Serious adverse events
    - Death
    - Acute myocardial infarction
    - Internal bleeding or hematoma
    - Arrhythmia including ventricular fibrillation
    - Low blood pressure/high blood pressure
    - Allergic reactions caused by contrast agent
    - Air/tissue/thrombotic embolus
    - Total occlusion of coronary artery
    - Allergic reaction to drugs
    - Coronary artery spasms
    - Bradycardia/palpitations
    - Catheter remnants left in body
    - Coronary artery occlusion/thrombosis

- Breakage, puncture, or injury to blood vessels
- Bleeding complications
- Peripheral occlusions
- Myocardial ischemia
- (Unstable) angina
- Infection and complications at puncture site
- Arteriovenous fistula
- Femoral pseudoaneurysm/pseudoaneurysm formation
- Cerebrovascular damage
- Other defects
  - Catheter kinking
  - Catheter twisting
  - Combined-use device insertion difficulty
  - Catheter cracking
  - Catheter insertion difficulty
  - Combined-use device removal difficulty

- Applications for pregnant, parturient and nursing women and children.

This device can be used in pregnant women or women who may be pregnant only if it is determined that the benefit of the treatment exceeds the risk of harm.  
[Because this device must be placed in the body under fluoroscopy.]

- Other cautions
 

Discard this product as medical waste after use.

#### Storage Instructions and Shelf-life

- Storage instructions
 

Avoid exposure to water, direct sunlight, and high temperature and humidity.
- Shelf-life
 

Refer to the expiry date indicated on the packaging.  
Shelf-life: 3 years after sterilization

#### Package

1 set / box

#### Symbols used for labeling

|  |                                  |
|--|----------------------------------|
|  | Manufacturer                     |
|  | Date of manufacture              |
|  | Use-by date                      |
|  | Batch code                       |
|  | Catalogue number                 |
|  | Sterilized using ethylene oxide  |
|  | Do not re-use                    |
|  | Consult instructions for use     |
|  | Non-pyrogenic                    |
|  | Do not use if package is damaged |

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