

HOCl
220 ppm

CLFine™

IONLESS™ HYPOCHLOROUS WATER

REVOLUTIONIZING WHOLE-ROOM DISINFECTION

BROAD-SPECTRUM | HYPOALLERGENIC | LOW-CORROSION | HIGH STABILITY



Hypochlorous Acid Water (HOCl)

BALANCING SAFETY, EFFICACY, MATERIAL COMPATIBILITY, AND SHELF STABILITY

Maintaining sterile environments while minimizing Healthcare-Associated Infections (HAIs) requires a disinfectant solution that harmonizes efficacy against pathogens, safety to humans, reduced corrosion, and shelf stability.

Traditional solutions like **bleach** - Sodium Hypochlorite (NaOCl) - may be shelf-stable and broadly effective, but they are also toxic to humans¹ and corrosive to medical equipment,² making them less than ideal for sensitive environments. The active ingredient in chlorine-based disinfectants like bleach is actually **hypochlorous acid (HOCl)**, a molecule produced naturally by the immune system, which is increasingly manufactured as a stand-alone disinfectant in many countries.⁴

HOCl is widely used as a topical antiseptic and wound-healing agent,¹⁴ as a disinfecting wash for vegetables and food processing equipment in schools and hospitals,¹⁶ and as a mouthwash for gargling at the dentist.²

Preferred by Healthcare Workers

HOCl has emerged as the superior chlorine-based alternative to bleach, due to its exceptional combination of safety, broad-spectrum efficacy, and ease of integration into primary and secondary disinfection protocols.^{2,7} Healthcare workers, especially in low-resourced hospitals, reported higher satisfaction and willingness to use HOCl, due to its odorlessness and reduced corrosiveness compared the strong odor and irritating fumes of bleach.^{5,6}


Profile Comparison of Hypochlorous Acid (HOCl) vs. Sodium Hypochlorite (Bleach)^{1,2,3,8}

	Hypochlorous Acid (HOCl)	Bleach (NaOCl)
Manufacturing	3 chamber electrolysis of diluted saline (NaCl) into HOCl and NaOH	Chemical reaction between chlorine (Cl ⁻) and sodium hydroxide (NaOH)
pH	5 - 6 (slightly acidic)	13 (highly alkaline)
Precautions	No special precautions necessary for atomization	Mask, rubber gloves, and goggles (with ventilation)
Instructions for Use	Programmable misting with CLeer 100	Bleach is a two-step process, must be rinsed with clean water to remove residue
Efficacy	HOCl reduces >99.9999% of bacteria ^{9,18} and renders >99.99% of viruses non-infectious in under 1 minute.	Bleach must remain wet on surfaces for up to 10 minutes for comparable efficacy to HOCl
Skin and Eye Safety	HOCl is gentle to the eyes and helps skin repair itself	Bleach can cause permanent nerve and tissue damage in the eyes and burns skin
Odor	Deodorizing	Strong
Side Effects	Hypoallergenic and non-irritating to the airways and mucous membranes	Bleach can cause nausea and severe irritation to mucous membranes
Disposal	Inactive HOCl degrades into a benign saline solution	Bleach requires hazardous waste disposal

Rapid Antimicrobial Action

Recognized by the World Health Organization (WHO) as a highly effective and hypoallergenic disinfectant, HOCl is safe for humans, environmentally friendly, and does not irritate the skin, eyes, or airways.¹⁶ HOCl's powerful, fast-acting antimicrobial action against bacteria, viruses, fungi, and even multidrug-resistant organisms is well-documented in clinical studies and laboratory tests.^{8,9,10,11,12,13}

Unlike negatively charged bleach which repels pathogens, HOCl has a neutral charge that attracts rather than repelling bacteria. Rather than slowly crushing pathogens by oxidizing their outer shell as bleach does, HOCl is absorbed by bacteria as if it were a nutrient, more easily penetrating biofilm¹³ and cell pathogen cell walls. This allows HOCl to exert an internal oxidizing effect directly on the DNA of viruses and bacteria^{8,9}, removing opportunities for resistance development.^{2,4,15}



Uncompromising Production Quality

While many HOCl products on the market are effective, they often come with compromises. Added sodium hydroxide (NaOH) buffers for stability can increase ion content, leading to greater conductivity and higher risks of surface oxidation and corrosion.^{2,3}

IONLESS™ Hypochlorous Acid

Nipro removes rather than adding impurities to stabilize its HOCl product, dramatically reducing its ion content and corrosivity. Purification through a reverse osmosis membrane produces an ultra-pure product having all the advantages of HOCl with none of the corrosion or stability issues seen with other everyday HOCl manufacturing brands.


Nipro produces a highly stable¹⁷ HOCl product with such extremely low ion content that not only is our product less corrosive than bleach, and less corrosive than other everyday HOCl products - **our HOCl product is less corrosive than tap water.**¹⁸

Ion Content and Conductivity of CLFine and Tap Water¹⁸

Contents	IONLESS™ HOCl	Tap Water in Japan	Tap Water Quality Standard
Effective chlorine concentration (ppm)	46	1	≤1
pH	6.14	7.73	5.8-8.6
Na ⁺ (ppm)	2.7	23.5	—
NaCl amount (ppm)	6.8	59.8	—
Anions (ppm)			
Cl ⁻	11.5	27.6	≤200
ClO ₃ ⁻	0.01	0.05	≤0.6
ClO ₂ ⁻	—	—	≤0.6
BrO ₃ ⁻	—	—	≤0.01
SO ₄ ²⁻	0.03	18.8	—
NO ₃ ⁻	0.03	4.67	≤10
Trihalomethane (ppm)			
CHCl ₃	0.005	0.015	≤0.06
CHBrCl ₂	0.0004	0.021	≤0.03
CHBr ₂ Cl	0.002	0.025	≤0.1
CHBr ₃	0.011	0.001	≤0.09
Total	0.018	0.061	≤0.1
Conductivity (μS/cm)	4.6	335	—

That is why we call it

IONLESS™



HOCl
220 ppm

Gentle to Humans, Deadly to Pathogens



IONLESS™ Hypochlorous Water

LESS CORROSIVE THAN TAP WATER & SHELF-STABLE,
DEACTIVATES THE MOST CONTAGIOUS AND DEADLY PATHOGENS

CLFine is an ultra-pure HOCl disinfectant that rapidly deactivates and reduces viral, fungal, and bacterial HAI²² with proven superior efficacy,^{19,20} while being hypoallergenic,² shelf-stable,¹⁷ and less corrosive than tap water.¹⁸

Paired with the CLeer™ 100 ultrasonic humidifier, CLFine 220ppm can disinfect entire unoccupied rooms at the push of a button, integrating easily into your hospital's disinfection protocols.

Broad-Spectrum

In suspension tests on the most common bacteria and viruses, HOCl at just 40 ppm concentration is equally or more effective than bleach at 1 000 ppm concentration in terms of efficacy and required contact time.²⁰

CLFine deactivates MRSA and Staphylococcus aureus, achieving a more than 99.9% reduction in adhered bacteria - at low concentration (40ppm) - 10 seconds in clinical suspension tests and within 4 hours in controlled atomization tests of 25 m² to 30 m² chambers.¹⁹

CLFine deactivates COVID-19 (SARS-CoV-2) rendering the virus non-infectious in under 20 seconds in clinical suspension tests from 30 ppm concentration and in under 60 minutes in controlled atomization trials.^{10,19}

Hypoallergenic, Humidifying, and Deodorizing

CLFine, unlike bleach, is gentle on human skin, non-toxic to the environment, and shows a high disinfecting and deodorizing effect.⁸

CLFine eliminates odors and provides safe humidification without the risk of harboring pathogens like Legionella.^{2,7}













Less Corrosive than Tap Water (IONLESS™)

Extended shelf-life achieved through reverse osmosis purification, ensuring prolonged effectiveness without compromising safety.

Accelerated material compatibility testing, under normal and humid conditions, found the corrosivity of CLFine 40 ppm equivalent to or less than the corrosivity of tap water.¹⁸

CLFine maintains a shelf life of 18 months when stored properly below 25°C and away from sunlight.¹⁷

Accelerated Testing for Material Compatibility¹⁸

Test material	CLFine IONLESS™ HOCl		Tap water	5% NaCl water
	airborne concentration 0.02 ppm	airborne concentration 0.5 ppm		
Stainless Steel	 Appearance change: None	 Appearance change: None	 Appearance change: Slight	 Appearance change: Slight
Iron	 Appearance change: Yes	 Appearance change: Yes	 Appearance change: Significant	 Appearance change: Significant
Electronic substrate component	 Appearance change: None	 Appearance change: None	 Appearance change: None	 Appearance change: Yes

Test sample exposure for 2 hrs (at 15-35°C). Allowed to stand in humid environment for 22 hrs (at 40°C/ 93%RH) x 4 days + Allowed to stand in standard atmospheric condition for 3 days (at 23°C/45-55%RH). (1 week per cycle x 5 cycles)

Antimicrobial Efficacy Rates of HOCl vs Bleach^{7,19}

Pathogens		Hypochlorous acid water (40 ppm)	Sodium hypochlorite (1 000 ppm)
Gram-positive bacteria	<i>Staphylococcus aureus</i>	⊙ (10 seconds)	⊙ (10 seconds)
	MRSA (methicillin-resistant staphylococcus)	⊙ (10 seconds)	⊙ (10 seconds)
	<i>Bacillus cereus</i>	△ (3-5 minutes)	△ (3-5 minutes)
	<i>Mycobacterium tuberculosis</i>	△ (2.5 minutes)	▲ (30 minutes)
	Other mycobacteria	△ (1-2.5 minutes)	▲ (2.5-30 minutes)
Gram-negative bacteria	<i>Salmonella Enteritidis</i>	⊙ (10 seconds)	⊙ (10 seconds)
	<i>Vibrio parahaemolyticus</i>	⊙ (10 seconds)	⊙ (10 seconds)
	<i>Escherichia coli</i>	⊙ (10 seconds)	⊙ (10 seconds)
	<i>Campylobacter jejuni</i>	⊙ (10 seconds)	⊙ (10 seconds)
	<i>Pseudomonas aeruginosa</i>	⊙ (10 seconds)	⊙ (10 seconds)
virus	<i>Candida albicans</i>	⊙ (10 seconds)	⊙ (10 seconds)
	Influenza viruses	⊙ (10 seconds)	⊙ (10 seconds)
	Coronavirus (SARS-Cov-2)	⊙ (10 seconds)	⊙ (10 seconds)
	Herpes virus	⊙ (10 seconds)	⊙ (10 seconds)
fungus	<i>Candida (Candida albicans)</i>	⊙ (10 seconds)	⊙ (10 seconds)
	Black mold (<i>Aspergillus niger</i>)	△ (5 minutes)	✗ (120 minutes)
	Blue mold (<i>Penicillium cyclopium</i>)	△ (5 minutes)	✗ (120 minutes)

Decreasing ratio of bacteria/viruses:

⊙ Not less than (NLT): 99.9% | ○ NLT 99% to LT 99.9% | △ NLT 90% to LT 99% | ▲ LT 90% | ✗ Ineffective





CLFine™
 IONLESS™ Hypochlorous Water

HOCl
 5 ppm

HOCl
 40 ppm

HOCl
 220 ppm

 **Cleer™ 100**

Ask us about Nipro's other CLFine products, including CLFine 5 ppm and 40 ppm, for misting in occupied spaces with Cleer™ 100.

- Automated sanitizing, humidifying, and deodorizing of waiting rooms, small event halls, and childcare playrooms.
- Cleans surfaces and clears away airborne pathogens with a gentle mist at low but effective concentration of CLFine 40 ppm.

 **Cleer™ 100**

NIPRO's Ultrasonic Humidifier for CLFine™

Intermittent Mode Setting

Setting	Mist	Stop
1	1 minute	3 minutes
2	1 minute	6 minutes
3	1 minute	9 minutes
User 1	0.5 to 90 minutes	0.5 to 90 minutes
User 2	0.5 to 90 minutes	0.5 to 90 minutes

User 1 and User 2 are customizable

Hybrid Mode Setting

1	2	3
10 minutes	20 minutes	30 minutes

Continuous mist timer auto-shifts back to intermittent



Continuous Mode
 3 levels of spray capacity



Speed Clean Mode
 Auto-shift to intermittent



Intermittent Mode
 Cycle of spray and interval



Key Lock / Off Timer
 For safety and efficiency

COMPLIANCE

Implementing Regulation (EU) 2021/347

EU regulation approves active chlorine released from hypochlorous acid (HOCl) as an active substance for use in biocidal products across several product types.

CDC Guidelines for Disinfection and Sterilization in Healthcare Facilities (2008)

CDC guidelines highlight hypochlorous acid as a highly effective and versatile disinfectant when used under appropriate conditions.

REFERENCES

1. Slaughter, Robin J et al. "The clinical toxicology of sodium hypochlorite." Clinical toxicology (Philadelphia, Pa.) vol. 57,5 (2019): 303-311. <https://doi.org/10.1080/15563650.2018.1543889>
2. World Health Organization. 2021. Hypochlorous Acid: Expert Review for Disinfection, Antisepsis, and Wound Care. WHO Expert Committee on the Selection and Use of Essential Medicines. 2021. Available at: https://cdn.who.int/media/docs/default-source/essential-medicines/2021-eml-expert-committee/expert-reviews/a18_hypochlorous-acid_rev2.pdf
3. Mehendale FV, et al. 2023. "HOCl vs. OCl⁻: clarification on chlorine-based disinfectants used within clinical settings." Journal of Global Health Reports. 2023;7:e2023052. <https://doi.org/10.29392/001c.84488>
4. da Cruz Nizer, Waleska Stephanie et al. "Surviving Reactive Chlorine Stress: Responses of Gram-Negative Bacteria to Hypochlorous Acid." Microorganisms vol. 8,8 1220. 11 Aug. 2020, doi:10.3390/microorganisms8081220
5. Gon, G. et al. 2022. A Better Disinfectant for Low-Resourced Hospitals? A Multi-Period Cluster Randomised Trial Comparing Hypochlorous Acid with Sodium Hypochlorite in Nigerian Hospitals: The EWASH Trial. Microorganisms 2022, 10, 910. <https://doi.org/10.3390/microorganisms10050910>
6. Naranjo-Soledad A et al. 2024. "Low-cost, local production of a safe and effective disinfectant for resource-constrained communities." PLOS Global Public Health 4(6): e0002213. <https://doi.org/10.1371/journal.pgph.0002213>
7. Functional Water Foundation. 2002. Guidelines for Hypochlorous Acid Water Generators (2nd Edition). <http://www.fwf.or.jp/kinousui.html>. Accessed 13 January 2024.
8. Severing, Anna-Lena et al. 2019. "Safety and efficacy profiles of different commercial sodium hypochlorite/hypochlorous acid solutions (NaClO/HClO): antimicrobial efficacy, cytotoxic impact and physicochemical parameters in vitro." The Journal of antimicrobial chemotherapy vol. 74,2 (2019): 365-372. <https://doi.org/10.1093/jac/dky432>
9. Hatanaka, N. et al. 2022. Hypochlorous acid solution is a potent antiviral agent against SARS-CoV-2, Journal of Applied Microbiology, Volume 132, Issue 2, 1 February 2022, Pages 1496–1502, <https://doi.org/10.1111/jam.15284>
10. Imoto, Y. et al. 2025. "Inactivation Effects of Hypochlorous Acid, Chlorine Dioxide, and Ozone on Airborne SARS-CoV-2 and Influenza A Virus." Food Environ Virol 17, 9 (2025). <https://doi.org/10.1007/s12560-024-09626-y>
11. Naka, A. et al. 2020. "Effectiveness of slightly acidic electrolyzed water on bacteria reduction: in vitro and spray evaluation." PeerJ. 2020 Feb 18;8:e8593. <https://doi.org/10.7717/peerj.8593>. PMID: 32110494; PMCID: PMC7034383.
12. Tsai, Cheng-Feng et al. 2024. "In vitro cytotoxicity and antibacterial activity of hypochlorous acid antimicrobial agent." Journal of Dental Sciences, Volume 19, Issue 1, 2024, Pages 345-356, ISSN 1991-7902. <https://doi.org/10.1016/j.jds.2023.07.007>
13. Palau M et al. 2022. "In Vitro and In Vivo Antimicrobial Activity of Hypochlorous Acid against Drug-Resistant and Biofilm-Producing Strains." Microbiol Spectr 10:e02365-22 <https://doi.org/10.1128/spectrum.02365-22>
14. Robson, Martin C et al. 2007. "Hypochlorous Acid as a Potential Wound Care Agent: Part II. Stabilized Hypochlorous Acid: Its Role in Decreasing Tissue Bacterial Bioburden and Overcoming the Inhibition of Infection on Wound Healing." Journal of Burns and Wounds vol. 6 e6. 11 Apr. 2007. PMID: 17492051
15. Dianty, Rahmi et al. 2023. "Electrolyzed hypochlorous acid water exhibits potent disinfectant activity against various viruses through irreversible protein aggregation." Frontiers in microbiology vol. 14 1284274. 19 Oct. 2023, <https://doi.org/10.3389/fmicb.2023.1284274>
16. Shen, L., Zhu, Y., & Chen, X. (2019). "Efficacy of Hypochlorous Acid in Reducing Microbial Contamination on Fresh-Cut Vegetables". Postharvest Biology and Technology, 152, 107-114. <https://doi.org/10.1016/j.postharvbio.2019.04.018>

Nipro Clinical Studies and Laboratory Tests

17. Internal Stability Testing
18. Internal Corrosion Study (Compatibility with Surfaces)
19. Internal Atomization Trials (Efficacy of Misting)
20. Internal Suspension Tests (Efficacy of Direct Application)

GLOBAL PRESENCE AND EXPERTISE

Nipro is a leading global healthcare company established in Japan in 1954. With over 38 000+ employees worldwide, Nipro serves the Medical Device, Pharmaceutical, and Pharma Packaging industries. As a total solution provider of renal care products, medical devices, and pharma packaging, Nipro supplies life-saving technologies, building long-term, meaningful partnerships with clinics and hospitals around the world. Together with its people, products, and processes, Nipro is committed to improving patient outcomes and quality of life, offering safe and superior products that optimize time, effort, and costs.

