Rebirth Pro2 Thrombus Aspiration Catheter

ENGLISH

WARNING

1. The product should only be used in facilities capable of performing emergency coronary artery bypass grafting (CABG) as a provision against complications that may cause injury and serious complications that could prove to be life threatening.
2. The product is a medical device and should only be used by physicians trained in percutaneous transluminal coronary angioplasty (PTCA).
3. Intravascular manipulation of the product should be carefully performed under fluoroscopy, confirming both tip movement and positioning. Movement of catheter tip or operation without confirming location may result in possible damage to vessel or the product.
4. Should any product resistance, abnormality in catheter tip movement or location occur during operation, immediately cease usage and confirm cause under fluoroscopy. Continued use may result in damage to vessel or the product.
5. When using other medical devices or pharmaceutical products together with the product, use them carefully and with a full understanding of their characteristics to prevent damage to the product.
6. Should any problem arise due to use in conjunction with another device, immediately cease use.
7. Care should be taken when undertaking drug administration through the product. Particles such as thrombus remaining in catheter may be flushed into patient.
8. Attention should be given to 'drop out' of thrombus during aspiration. 'Drop out' refers to the phenomenon when aspirated thrombus becomes caught at catheter tip, guide catheter tip or at stenotic lesion site. Performing angiography under this condition forces the thrombus into the distal vasculature. Should drop out be suspected, remove blood contained within the guiding catheter. For other methods to overcome this situation, please refer to section titled PRECAUTIONS.
9. Reference package insert for each medical device and pharmaceutical product used together with this product.
10. During procedure, perform administration of the appropriate anticoagulants and vasodilators, as directed by lead physician.
11. All product use should be undertaken in a sterile environment.

CONTRAINDICATIONS

1. Contraindications for Use

- (1) The product is sterilized, is not re-usable and must not be re-sterilized. Resterilization and/or reuse could result in infection or degradation of product characteristics such as shaft strength or lubricity and could result in failure of the product during use.
- (2) This product is only for intravascular use.

2. Lesions

Lesions in left main trunk that is not protected by collateral blood flow or bypass.

3. Patients

- (1) Not appropriate for Coronary bypass surgery.


- (2) Unable to accept appropriate sheath size due to anatomical abnormality.
- (3) Exceedingly long coagulation time or unable to undertake anti-platelet, anti-coagulation therapy.
- (4) Severe allergic reaction to any pharmacological agent used during the course of a normal PCI procedure.
- (5) Low left ventricular function.
- (6) Bacteremia or blood poisoning.
- (7) Hemorrhagic complications.
- (8) Pregnant or suspected of being pregnant.
- (9) Unstable hemodynamics or shock.
- (10) Severe disease of other internal organs.
- (11) Determination by physician that use of the product is not appropriate for reasons other than those noted above.

4. Prohibition of Use with Medicine, Medical Equipments

- (1) Do not use any products comprising of alcohol, organic solvent, fat-based emulsions or oil to prevent damage to the product and/or hydrophilic coating deterioration.
- (2) Do not undertake contrast media administration through product from an injector system. The injection pressure may cause damage to the product or cause residual thrombus within the product to be flushed into the vasculature.

SHAPE / CONSTRUCTION

1. Product Diagram

2. Accessories Diagram


The product packaging includes product and accessories. The contents of accessories include: flushing device, syringe, stopcock (two way), stopcock (three way), extension tube, cell strainer, filter, drainage tray and stylet.

- (1) Stylet


NOTE: The stylet has been pre-inserted in the product (aspiration catheter)


- (2) Syringe


- (3) Extension tube

- (4) Stopcock (two way)

- (5) Stopcock (three way)


- (6) Flushing device

- (7) Cell strainer

- (8) Filter and drainage tray

PURPOSE, EFFICACY

- For percutaneous removal or dissipation of thrombus from the peripheral vasculature.
- For aspiration of coronary thrombus.

SPECIFICATIONS

- Tensile Strength Of Catheter Junction:
 - 40 mm section between Proximal and Distal Shaft: $\geq 5N$
 - 40 mm section between Hub and Proximal Shaft: $\geq 5N$
- Conformity Of Connector: ISO594-1, 594-2

METHOD OF USE

1. Preparation prior to Use

- The applicable sizing of the guiding catheter and the guide wire to be used in conjunction with the product are as noted below:

Rebirth Nominal Size	Applicable ID of Guiding Catheter (minimum diameter)	Applicable Guide Wire (maximum diameter)
6Fr	0.070" (1.78mm)	0.014" (0.36mm)
7Fr	0.080" (2.03mm)	0.014" (0.36mm)

ID : Inner Diameter

- Inspect and confirm each device for correct functionality.
- Remove the product from the packaging (still in the protective hoop).
- Flush the product with heparinized saline while retaining the product in the protective hoop. Ensure the product does not become contaminated by being ejected from the hoop during flushing process.
- Carefully remove the product from the protective hoop.
- Connect the extension tube to the connector fitting on the stylet, and the stopcock to the extension tube. When either or both extension tube and stylet are not required, directly connect the stopcock to the product or the stylet connector.
- Insert 10ml of heparinized saline into the syringe and remove the air contained within the syringe. As per Fig.1,

attach the stopcock to the syringe and open the stopcock valve.


Fig.1:

Example of Combination / Connection of Accessories

- Using the saline contained within the syringe, flush the entire amount through the catheter aspiration lumen including the stopcock (and the extension tube when applicable), repeating the process should residual air still remain and closing the stopcock when process is completed.
- Attach the flushing device to a separate syringe containing saline and completely flush the guide wire lumen through from the product tip (tip of the guide wire lumen).

2. Insertion of Product

- In accordance with the Instruction For Use of any pharmaceuticals or devices used in combination with the product, undertake preparation of the product prior to insertion.
- Place the end of guide wire into the product tip and insert the product until the product tip meets the guide wire tip.
- Carefully insert the product through the hemostatic valve connector and while confirming under fluoroscopy advance the product until the radiopaque marker reaches the target location.
- When closing the valve of the hemostatic valve connector to minimize blood loss, do not over tighten for smooth product insertion and avoiding damage to the product.
- Remove the stylet from the product (if the stylet has been utilized) disconnecting the extension tube from the stylet, and reconnect the extension tube to the product.

3. Thrombus Aspiration

- After confirming location of the product under fluoroscopy, pull the plunger of the syringe back to the required volume and lock it in position under negative pressure.
- Opening the stopcock, slowly advance and retract the product at location of thrombus while aspirating.
- After thrombus aspiration, return the product to the proximal location to where thrombus existed.
- After closing the stopcock, remove the syringe.
- Should further aspiration be required, remove the aspirated material contained within the syringe, re-connect to the product and repeat steps (1) to (3).

4. Removal of Product

- Confirm completion of thrombus aspiration.
- Confirm the stopcock is in a closed position [removal of product without closing the stopcock may cause the introduction of air into the guiding catheter].
- Slowly remove the product.
- Confirm there is no residual thrombus within the hemostatic valve connector. Remove the remaining thrombus from the hemostatic valve connector.
- If the product is to be reintroduced into the patient, re-insert stylet (as required) within the product and commence procedure from re-flushing.

PRECAUTIONS

1. Precautions prior to Use

- (1) Confirm specifications of the product and other devices used in conjunction for compatibility.
- (2) Do not use the product if the packaging or contents are damaged or contaminated.

2. Precautions during Use

- (1) Always precede insertion of the product with initial placement of the guide wire. Inserting the product without the guide wire may potentially cause damage to arterial wall or dissection.
- (2) Ensure that the guide wire has been placed adequately distal to the target lesion. Failing to place the guide wire in a distal location may cause the guide wire to drop out of the guide wire lumen, potentially causing damage to vessel.
- (3) Should guide wire drop out of guide wire lumen of product, remove entire product from patient and re-insert guide wire.
- (4) There may be product resistance should the guide wire become entangled. Continuing to insert the product without correcting entanglement may result in damage or breakage to either the product and/or the guide wire.
- (5) Slack in the guide wire may cause a reduction in handling performance. The guide wire may potentially become stuck at the tip of the guiding catheter when retracting the product given the guide wire becomes loose at the guide wire exit port location, especially for instances where slack occurs between the guiding catheter and guide wire exit port of the product. Remove excess slack by slowly moving the guide wire forward and backward. Damage or breakage may occur if the product or the guide wire is handled forcefully.
- (6) When crossing a previously implanted stent with the product, advance the product ensuring that the product tip does not become caught on the stent. Damage to either the stent or the product may result.
- (7) Care should be taken not to cause damage to the product by other utensils (scalpel, scissor etc). Should damage occur, immediately cease use of the product.
- (8) Remove all air particles contained within the product and the syringe. Adequate aspiration may not be achieved given residual air.
- (9) The product should only be used after being adequately soaked in heparinized saline. Failure to adequately flush the product may result in the coagulation of blood or contrast media within the guiding catheter, as well as reduce the lubricity of the hydrophilic coating and negatively affect handling.
- (10) Do not force insertion of the product or retrieve the product rapidly for highly tortuous sections of artery, bifurcations or calcified lesions. The product may kink or be damaged, as well as result in damage to the artery.
- (11) Operation of the product should be performed while confirming location of thrombus under fluoroscopy. Advancing the product distal to location of thrombus may cause thrombus to flow to the distal bed, and the product tip may aspirate the arterial wall during aspiration causing damage to the arterial wall.
- (12) The product should be inserted and removed carefully. Rapidly moving the product may cause kinking of the product, damage to the guide wire lumen or cause thrombus to flow into the distal bed.
- (13) Should the product become bent, broken or kinked, immediately cease use. Proceeding to use the product in

this condition may decrease aspiration function, further damage the product or potentially leave separated parts of the product within the artery.

- (14) Do not apply torque to product. Rotating the product may cause damage or breakage of the product, potentially leaving residual components within the vessel.
- (15) Place the product super-selectively by deep seating the guide catheter. This avoids unnecessary damage to arterial wall when inserting the product, as well as reduce the potential for thrombus within the guiding catheter being pushed into the peripheral vasculature.
- (16) Ensure that the stylet tip does not damage the product when reinserting the stylet.
- (17) Ensure that aspirated blood, thrombus or air contained within the syringe is not administered from the syringe through the product.
- (18) Ensure that aspirated blood is not re-administered into the patient. Infection or possible adverse events may result.
- (19) When placing the syringe under negative pressure, take care not to catch fingers between the plunger and the syringe barrel.
- (20) Appropriate volume of aspirated thrombus is determined by acting physician, with the appropriate volume made by adjustments to syringe.
- (21) Do not advance the product without undertaking thrombus aspiration. There is a potential for thrombus pushed by the product to travel into the distal vasculature causing embolization.
- (22) Always undertake thrombus aspiration with the stylet removed from the product to avoid a blockage with the thrombus causing ineffectual aspiration.
- (23) Ensure aspiration is performed to avoid dissipated thrombus from traveling into the peripheral vasculature or side branch. Sudden operation may cause thrombus to embolize peripheral vasculature.
- (24) Should interruption of flow be confirmed during thrombus aspiration, slowly pull back the product until flow has been reestablished, and re-advance the product tip. Should flow not be reestablished, there is potential for thrombus to be stuck in the product or at the product tip. Remove the entire system under the condition of negative pressure and flush the entire system. Positive pressure through the product while in-vivo may potentially cause thrombus to flow to distal vasculature.
- (25) Should interruption of flow be confirmed during thrombus aspiration, confirm that it is not the arterial wall being aspirated. Further, confirm that the guiding catheter is not in a wedged condition and is running parallel to the coronary artery.
- (26) Should interruption of flow be confirmed during thrombus aspiration, completely remove the entire system under a condition of negative pressure. Removing only the product may cause thrombus to drop out.
- (27) When carefully retrieving the product through the guiding catheter with blockage of product tip, take care to ensure that such thrombus does not drop out. When it is suspected that the thrombus dropped out, be sure to retrieve blood contained within the guiding catheter prior to undertaking further actions to avoid vascular embolization.
- (28) When removing the product (excluding when thrombus is caught within the guiding catheter), always ensure that the stopcock is in a closed position. If the stopcock is not in a closed position, there may be potential that the aspirated thrombus may be re-introduced into vessel,

or blood within the guiding catheter may be aspirated, potentially causing an introduction of air into the guiding catheter.

- (29) During removal of the product, fix the proximal section of the guide wire with hand. The guide wire may be potentially removed together with the product if not fixated. Further, trying to reintroduce the product without fixating the guide wire will result in dropping out the guide wire from the product.
- (30) Do not flush contrast media or introduce catheters through the hemostatic valve connector connected to the product should thrombus remain in the hemostatic valve connector, because thrombus may potentially flow to vessel.

3. Precautions after Use

Dispose of this product as medical waste, and take measures to prevent possible spread of infection.

4. Adverse Event

Possible adverse events include, but are not limited to, the following:

Death, myocardial infarction, internal hemorrhage, hematoma, ventricular fibrillation including arrhythmia, hypertension, hypotension, hemorrhagic complications, arterial spasm, stroke, distal embolization, arterial or bypass graft occlusion, arterial rupture or perforation or damage, unstable angina, medicinal reaction or allergic reaction to contrast media, infection, arteriovenous fistula, air embolization, arterial dissection, blood loss from puncture site, ischemia, caused by long duration inflation, intravascular thrombosis, nausea or vomiting, palpitation, tachycardia, bradycardia.

STORAGE METHOD, SHELF LIFE AND OTHER

1. Storage Method

- (1) Store the product in a room temperature location not exposed to high temperature and humidity or direct sunlight and take proper precautions to ensure product does not contact water.
- (2) Avoid inclinations, vibrations and impacts (including during transportation) and store in a safe, stable environment.
- (3) Do not store near chemicals or in areas where the device may be exposed to gases.

2. Shelf Life

Use this product before the "expiration date" shown on the package label.

3. Sterilization Method

The product package has been sterilized by exposure to Ethylene Oxide Gas (EtO) and not intended to be re-sterilized.

PRODUCT WARRANTY DISCLAIMERS AND LIMITATION OF REMEDIES

WITH RESPECT TO ANY AND ALL GOODMAN PRODUCT(S) DEPICTED OR DESCRIBED IN THIS PUBLICATION, THERE ARE NO EXPRESS OR IMPLIED WARRANTIES OF ANY KIND WHATSOEVER, INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. GOODMAN CO., LTD AND ITS SUBSIDIARIES ("GOODMAN") SHALL NOT BE LIABLE FOR ANY DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES IN CONNECTION WITH THE USE,

RE-USE, OR ANY OTHER ASPECT OF THE RELEVANT PRODUCT(S) EXCEPT AS EXPRESSLY PROVIDED UNDER APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND GOODMAN TO ANY REPRESENTATION OR WARRANTY. ANY PRODUCT OR SERVICE DESCRIPTION OR SPECIFICATIONS CONTAINED IN ANY GOODMAN PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE RELEVANT PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS OR OTHER WARRANTIES.

	Catalogue number
	Batch code
	Use by
	Date of manufacture
	Do not reuse
	Do not resterilize
	See instructions for use
	Sterilized using ethylene oxide
	Keep dry
	Keep away from sunlight
	Do not use if package is damaged
	Manufacturer

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Goodman Co., Ltd.