NSE PTA Balloon Dilatation Catheter

EN	English	English	2
DE	German	Deutsch	8
FR	French	Français	11
ES	Spanish	Español	14
IT	Italian	Italiano	17
SV	Swedish	Svenska	20
FI	Finnish	Suomi	23
NL	Dutch	Nederlands	26
PT	Portuguese	Português	29
CS	Czech	Česky	31
SK	Slovak	Slovensky	35
PL	Polish	Polski	38
HU	Hungarian	Magyar	41
ZH-TW	Chinese (Traditional Chinese)	中文 (繁體中文)	44
RU	Russian	Русский язык	46
	Labeling Symbol Definitions		50

MANUFACTURER:

Goodman Co., Ltd.
 5F KDX Nagoya Sakae Building,
 4-5-3 Sakae, Naka-ku, Nagoya,
 Aichi 460-0008
 Japan

ENGLISH

***This English IFU applies for countries other than Japan and China.**

WARNING

1. Safety and efficacy has not been clearly determined for the use of the device when applied to complex and tortuous lesions in the iliac artery. Careful consideration should be made in regards to the potential for vessel rupture requiring surgical treatment and risk of occlusion of the distal artery from vasculature dissection.
2. Do not exceed balloon expansion beyond the vessel size both proximal and distal to the lesion (damage to vessel may result).
3. As the elements are fixed both ends, care should be taken to ensure that additional devices are not passed through the gap between element and balloon (the secondary device may become trapped, damage to either product may result).
4. Care should be taken when inserting and expanding device at calcified and in-stent restenotic lesions and at location of stiff synthetic grafts (damage to vessel or product may result).
5. Care should be taken to ensure that the device does not become entangled with multiple stents implanted at bifurcated lesions (damage to product may result).
6. Care should be taken when used at the distal location to a drug eluting stent (disruption of newly formed neointimal formation may result).

7. The device should only be manipulated in a completely deflated position (damage to product may result).

CONTRAINDICATIONS

1. Contraindications For Use


The product is sterilized, is not re-usable (single use only) and must not be re-sterilized. Resterilization and/or reuse could result in infection or degradation of product characteristics such as balloon size, shaft strength or lubricity and could result in failure of the product during use. This device is not intended for use for post dilatation of a stent.

2. Inappropriate Lesions

- (1) Exceeding the stent strut in a bifurcated lesion
- (2) Within a newly implanted stent or distally to the stent
- (3) Fractured stent
- (4) Lesions in coronary and intracranial vessels including carotid arteries

SHAPE / CONSTRUCTION

1. Product Diagram

* Devices with a working length of 50cm and 90cm do not have depth markers.

<Components>

Catheter clip (Devices with a working length of 50cm do not contain the catheter clip.) / 1 unit

<Material>

Nylon resin, polyethylene resin, polycarbonate, polyethylene oxide.

2. Product Summary

The product is an OTW PTA balloon catheter used for the purpose of dilatation of stenotic lesions during percutaneous transluminal angioplasty.

The balloon expands to the specific OD and length at nominal pressure and has three elements on the outside of balloon (with elements creating a focused transmission of dilative force, i.e. scoring effect) in order to dilate stenosis that is traditionally considered difficult to expand.

This product has a hydrophilic coating applied to the surface.

PURPOSE, EFFICACY

This device is used for percutaneous transluminal angioplasty (PTA) to dilate a stenotic lesion in peripheral arteries and stenotic lesions of arteriovenous dialysis fistula, whereby it is considered that dilatation of a conventional balloon would result in inadequate lesion expansion. This device is not for use in the coronary and intracranial vessels including carotid arteries.

APPROPRIATE PATIENT SELECTION

Patients with appropriate peripheral arteries for a 4-7mm balloon diameter in a percutaneous transluminal angioplasty procedure as determined by a physician (including pregnant patients).

SPECIFICATIONS

1. Tensile Strength Of Connections

- 1) Proximal shaft – hub: $\geq 10.0N$
- 2) Balloon – distal shaft: $\geq 10.0N$
- 3) Balloon – distal tip: $\geq 3.0N$
- 4) Cover tube – distal tip: $\geq 3.0N$

2. Nominal Pressure (NP)

8atm [$8 \times 10^2 kPa$]

3. Rated Burst Pressure (RBP)

- 1) Balloon diameter 4.0mm:
18atm [$18 \times 10^2 kPa$]
- 2) Balloon diameter 5.0, 6.0mm:
16atm [$16 \times 10^2 kPa$]
- 3) Balloon diameter 7.0mm:
14atm [$14 \times 10^2 kPa$]

*In vitro test results indicate that 99.9% of the balloons do not incur rupture (95% reliability) when expanded up to maximum recommended inflation pressure.

4. Conformity Of Connector

ISO80369-7

MEDICAL DEVICES REQUIRED FOR THIS DEVICE

- Guidewire
Maximum guide wire diameter:
0.46mm (0.018inch)
- Sheath introducer
Minimum Introducer Sheath
 - 1) Balloon diameter 4.0, 5.0mm: 5F
 - 2) Balloon diameter 6.0, 7.0mm: 6F
- Inflation device with manometer
- Stopcock with luer connection conform to ISO80369-7
- Luer-lock syringe
- Hemostatic valve connector

METHOD OF USE

1. Preparation

- (1) Inspect and confirm each device for correct functionality.
- (2) Remove the product from the packaging (remaining in the protective hoop).
- (3) Carefully remove the product from the protective hoop.
- (4) Carefully remove the balloon protector and the stylet. Confirm that there is no damage to the product.

*NOTE: Inflation and deflation of the balloon shall not be performed prior to use unlike the preparation for a conventional PTA balloon catheter.

- (5) Ensure that entire length of shaft of device is saturated in heparinized saline solution.
- (6) Connect the stopcock to the product hub.
- (7) Prime an inflation device with appropriate amount of diluted contrast media (contrast: heparinized saline = 1:1). Connect it to the stopcock and flush the port of stopcock with contrast media.
- (8) With the tip of the inflation device pointing downwards, remove all residual air within the balloon lumen by continuous application of negative pressure.
- (9) With negative pressure applied to the inflation device, close the stopcock ensuring that negative pressure of the balloon lumen is maintained.

*Do not allow air or liquid to enter the balloon in order to preserve the balloon folding functionality and the elements.

- (10) Attach the syringe to the hub guide wire port, flush the guide wire lumen of the device with heparinized saline, removing all residual air.

2. Insertion Of Product

- (1) Follow the package insert for each medical device and complete preparation of the devices prior to insertion of the product.
- (2) With the balloon completely deflated, insert the guide wire into the distal tip of the product and carefully advance the product to the target lesion.

3. Balloon Inflation

- (1) Confirm position of the balloon in relation to the target lesion and close the hemostatic valve, locking balloon into place.
- (2) While confirming the balloon diameter under fluoroscopy, inflate and deflate the balloon using an inflation device (refer to the compliance chart included in packaging for each balloon size reference diameter at corresponding inflation pressure).
- (3) Perform multiple inflations as deemed necessary.
- (4) Upon completion of dilatation, ensure the balloon is completely deflated, return the product to within the guiding catheter and evaluate improvement to stenotic site via angiography.

4. Removal Of Product

While confirming guide wire positioning, carefully remove the device through the sheath introducer (guiding catheter).

PRECAUTIONS

<Fundamental Considerations>

- (1) This device should be used in those facilities capable of immediate surgical response to any complication with potential negative effects on health or life threatening complication that may arise as a result of its use.
 - (2) Use of this product should be undertaken under an appropriate anticoagulation regimen.
 - (3) Appropriate balloon sizing and products used in combination should be determined based on the anatomical structure of the patient. Size selection (both diameter and length) should not exceed the proximal and distal sites of the stenotic location, nor exceed lesion length.
 - (4) All device manipulation should be carefully performed under angiography while constantly confirming the position of the catheter (damage to vessel and product may result).
 - (5) Should any resistance be incurred during usage or any anomaly of the catheter position be observed, immediately cease product manipulation and confirm the cause under angiography (damage to vessel or product may result).
 - (6) Should any damage to device be observed, immediately cease use and confirm condition of the elements (damage to vessel may result).
 - (7) All device removal from protective packaging and hoop, as well as removal of sheath and stylet should be performed carefully so as to not place excessive burden on the product (product insertion and functionality may be affected).
 - (8) Use a gauze that has been sufficiently dampened with heparinized saline to remove any particles attached to the guide wire when inserting and removing the device (operability of guide wire and product may be affected, damage may result).
 - (9) This device should only be used by physicians trained in percutaneous intravascular procedures.
 - (10) Do not exceed RBP (damage to device may result).
 - (11) All balloon inflations and deflations should be performed under radiocopy while confirming the entire balloon.
 - (12) Do not rotate device (damage may result).
 - (13) Do not use products comprising of organic solvent, fat-based emulsions or oil (damage to hub and/or elements, as well as deterioration to the hydrophilic coating may result).
 - (14) Careful consideration of pharmacological allergies to products used, potential of patient for being pregnant, hemodynamic condition and potential for cardiogenic shock should be undertaken.
-
- 1. Precautions Prior To Use**
 - (1) Use of distal protection is recommended for those patients with potential for distal embolization.
 - (2) Refer to attached documents of all pharmaceutical and medical devices used in combination.
 - (3) Do not use the product if the packaging or contents are damaged or contaminated, in particular that dam-

age to elements has not incurred.

- (4) All product use should be undertaken in a sterile environment.

2. Precautions During Use

- (1) Do not cause damage to product when using sharp items.
- (2) All product manipulation should be performed with the guide wire advanced from device tip.
- (3) All balloon inflations should be performed with contrast media, ensuring there is no residual air (potential for air embolization may result).
- (4) During balloon inflation confirm condition of inflation device as well as inflation via fluoroscopy. Should any irregularity be observed, immediately cease inflation, deflate balloon and remove the device, confirming potential causes.
- (5) During long periods of use carefully consider thrombus accumulation (device and guide wire functionality may be affected).
- (6) During inflation in a position that does not allow for concentric balloon expansion, care should be taken to ensure the balloon does not move position (damage to vessel may incur).
- (7) Do not overtighten the hemostatic valve, potentially affecting guide wire manipulation, as well as balloon inflation and deflation by impeding flow of contrast media.
- (8) Do not use product if kinking is observed (product may break).
- (9) When using multiple devices in-vivo, ensure that devices do not become entangled. Should any resistance during device delivery be encountered,

immediately determine the cause (product damage may result).

- (10) Continuously confirm the condition of the catheter to ensure that there is no damage, that connections are attached and that no leakage of contrast media is determined. Should any abnormality be determined immediate cease use of the product (complications to the vessel may result).
- (11) Once the device has been removed from the balloon protective sheath, do not re-insert the device (damage to product may result).
- (12) This device has a hydrophilic coating and must remain damp from heparinized saline at all times (operability may be affected and damage may result).

3. Precautions After Use

Dispose of the product as medical waste and take measures to prevent possible spread of infection.

4. Adverse Events

Possible adverse events include, but are not limited to, the following:

- death
- myocardial infarction
- acute closure
- closure, ischemia
- arrhythmia including ventricular fibrillation
- angina
- cerebrovascular disease
- restenosis
- claudication
- necrosis
- quadruple amputation
- embolization (air, tissue, thrombus)

- complications of hemorrhage
- puncture complications
- hematoma
- pain
- vascular spasm
- vessel dissection, perforation, rupture
- arteriovenous fistula
- pseudoaneurysm
- infection
- pharmaceutical allergic reaction

STORAGE METHOD, SHELF LIFE AND OTHER

1. Storage Method

- (1) Store the product in a room temperature location not exposed to high temperature and humidity or direct sunlight and take proper precautions to ensure product does not contact water.
- (2) Avoid inclinations, vibrations and impacts (including during transportation) and store in a safe, stable environment.
- (3) Do not store near chemicals or in areas where the device may be exposed to ionizing radiations.

2. Shelf Life

Use this product before the date of "Use By" shown on the package label.

3. Sterilization Method

The product is supplied STERILE using an Ethylene Oxide Gas(EtO).

GOODMAN PRODUCT(S) DEPICTED OR DESCRIBED IN THIS PUBLICATION, THERE ARE NO EXPRESS OR IMPLIED WARRANTIES OF ANY KIND WHATSOEVER, INCLUDING, AND WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. GOODMAN CO., LTD AND ITS SUBSIDIARIES ("GOODMAN") SHALL NOT BE LIABLE FOR ANY DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES IN CONNECTION WITH THE USE, RE-USE, OR ANY OTHER ASPECT OF THE RELEVANT PRODUCT(S) EXCEPT AS EXPRESSLY PROVIDED UNDER APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND GOODMAN TO ANY REPRESENTATION OR WARRANTY. ANY PRODUCT OR SERVICE DESCRIPTION OR SPECIFICATIONS CONTAINED IN ANY GOODMAN PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE RELEVANT PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS OR OTHER WARRANTIES.

PACKAGING

1unit / box

PRODUCT WARRANTY DISCLAIMERS AND LIMITATION OF REMEDIES

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