

# GravineXt™


ACID CONCENTRATE DELIVERY SYSTEM



# GravineXt™

## ACID CONCENTRATE DELIVERY SYSTEM

Designed to transfer dialysis concentrate from storage tanks to the dialysis machine.



### Working Principle

1. The pump moves the concentrate from its tank to the gravity tank.
2. From the gravity tank, the concentrate is delivered to dialysis machines by gravitational force.
3. Sensors in the gravity tank automatically detect when it must be refilled.

### Flexible Installation

4. GravineXt can be connected to up to 4 tanks, with up to 2 different types of concentrate. For each type of concentrate, while one tank is in operation the other is on standby.

### Versatile Supply

5. GravineXt is suitable for tanks, octabins (with top or bottom connections), bag-in-boxes, mixing systems of any sizes, and all types of formulas.

### Monitoring GravineXt

6. Control panel with 5.7" (14 cm) touch screen:
  - in the storage room
  - on the dialysis unit (optional)
  - displayed on a computer and/or phone (Ethernet/SIM card)

# Why use gravitational force?


Gravity ensures equal pressure across the loop, reducing the risk of leakage.



## Ethernet Connection

- Warning when octabin/bag-in-box is empty
- Technical service team notified of failures
- Actions and alarm failure history

# Why choose GravineXt?



## Uninterrupted supply in automatic mode

- Automated switching between concentrate tanks
- Optional automatic start with reverse osmosis systems
- Level indicator for volume remaining in the octabin/bag-in-box

## Quality & Safety

- Double pump system per concentrate type, allowing high operational safety in the unlikely event of one pump failing
- Dry-running safety pumps
- Uses corrosion-free polypropylene
- Low pressure system
- Four level sensors
- Alarm log

## Simple & Flexible


- Single-use concentrates can be completely replaced
- Reduces dialysis preparation time
- Minimal ecological footprint
- Frees up storage space

# Technical Data

## Performance

Max. pump capacity	5.7 L/min
Max. height from pump module to gravity tank	30 m
Max. distance from pump module to gravity tank	50 m
Gravity tank capacity	54 L or 80 L
Max. beds	30 (54 L)   60 (80 L)

Max. number of beds



## Electrical data

Supply voltage	230 V / 50 Hz (± 10%)
Max. electrical consumption	0.9 kW
Pre-fuse	16 A
	Leakage current protection 30 mA

## Connections

Supply tanks connection	Hose connection DN 10x8, connectors
Loop (gravity tank to dialysis machines)	Hose connection DN 10x14, connectors
Electrical connection	CEE 7 standard

## Ambient temperature

Storage / transport	1-35°C
Operation	10-35°C
Relative humidity	< 95% without condensation
Air pressure	795-1062 hPa

## Dimensions


Switch cabinet	480 x 690 x 240 mm
Gravity tank (54 L)	700 x (200 + 55 (vent filter)) x 400 mm
Gravity tank (80 L)	700 x (300 + 55 (vent filter)) x 400 mm
Pump module	710 x 400 x 160 mm



## Make it easy

Move from individual containers to an automated central supply with GravineXt

### Save on effort



handled by nurse

## Make it greener<sup>1</sup>



**92%\***

reduction in water use from the production of acid and its components

**Save on water**



**40%\***

Reduction in kg of CO<sub>2</sub> eq<sup>2</sup>

**Save on emissions**



**48%\***

reduction in carton waste

**Save on packaging**



**78%\***

reduction in volume of plastic waste

## in pursuit of a better tomorrow

1. French analysis of GravineXt life cycle compared to single-use bags

2. Evaluation by In Extenso organization. 2022. Report available upon request.

\* up to



**NIPRO**  
Live Longer. Live Better.

**Nipro Medical Europe**

Blokhuisstraat 42, 2800 Mechelen, Belgium

[www.nipro-group.com](http://www.nipro-group.com) | [medical@nipro-group.com](mailto:medical@nipro-group.com)