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en: Instructions for use Acid Concentrate (acetate based) and Citrasate™ (citrate based) for Bicarbonate Haemodialysis
Indications: Preparation of dialysis fluid for the treatment of patients with acute or chronic renal failure when haemodialysis or related therapies are prescribed by a physician. **Intended use:** Liquid acid concentrate for bicarbonate haemodialysis. Use only with bicarbonate haemodialysis concentrate in liquid (NaHCO₃ 8.4% or 6.6%) or powder form in the given dilution. If the haemodialysis monitor is intended to be used for online treatments, the liquid acid concentrate can also be used in this setting (see relevant machine operator's manual). **Side effects:** No specific side effects have been reported for the acid concentrate. However, the dialysis treatment can have the following side effects: hypotension, hypertension, arrhythmias, nausea, vomiting, spasms, chills, headache, dizziness, numbness, palpitations, fatigue, malaise, itching, alkalosis, reduced cerebral blood flow, respiratory distress, and inflammation have been reported. Bleeding or clotting during treatment may result in low dialysis efficacy, haemorrhage, and vascular access problems. **Contraindications:** Patients who are hypersensitive to acid concentrates (A-component (Acetic acid) and Citrasate™ (Citric acid)) or any of their components should not be dialyzed with this product. The type of acid concentrate formula used is dependent on the patient's electrolytical status. The physician must select the correct formula based on the patient's requirements, to avoid for

MTN Neubrandenburg GmbH
Gustav-Kirchhoff-Str. 2
17033 Neubrandenburg, Germany
www.nipro-group.com

2797

0086

UK Responsible Person:
Nipro Medical UK Ltd.
Units 12-14 South Point,
Ensign Way, Hamble,
Hampshire
SO31 4RF, United Kingdom

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example hyperkalaemia (applicable for potassium-containing concentrates), hypokalaemia (applicable for potassium-free concentrates), inadequate magnesium or calcium-plasma concentration, and unmanageable coagulation disorders. For patients with strong haemodynamic instability, an alternative method of extracorporeal treatment may be indicated based on the overall patient indications. **Warnings:** Overmixing of bicarbonate concentrate must be avoided. The concentration of the ready-to-use dialysis fluid must be precisely monitored. Before dialysis treatment and each time a new concentrate is used with the machine, check the conductivity and pH-level of the dialysis fluid. Acid and bicarbonate concentrate should only be diluted with water of suitable quality according to ISO 23500-3 immediately before use. Recommendation for the use of opened containers: small containers (canisters, bags, poltainer) max. 1 week, large containers (Container, Bag in Box) max. 1 month. Do not use large containers (Container, Bag in Box) after disconnection from central ring line system. Do not use after the expiration date. Only use undamaged containers and clear colorless liquid. **Note:** Bag for use with Nipro haemodialysis monitors only. Follow the dialysis machine manufacturer's instructions. Depending on the volume, the concentrate can be used for more than one dialysis treatment. Keep the container closed when not in use. For safe connection of Containers or Bag in Box to central ring line systems, the instructions in the appendix to the instructions for use must be followed. Concentrate is endotoxin-free. For safe application and disposal, see the manufacturer's safety data sheet. Any serious incident that has occurred in connection with the product must be reported to the manufacturer and the competent authority of the state in which the user is established. **Notes on Citrasate™:** For better differentiation from acid concentrate with acetic acid, the Citrasate™ label is marked with a yellow bar. Check metabolizing capacity of citrate for all patients. Patients with low muscle mass and/or restricted liver function may not be suitable without further correction of concentration of free calcium concentration. Compared with acetate dialysate, notably higher calcium and magnesium concentrations may be required in the citrate-containing dialysate to achieve similar systemic post-dialytic values. Because one citrate ion will be converted quickly during dialysis treatment into 3 bicarbonate ions, a change of acid/base balance occurs with the possible result of temporary alkalosis. Therefore, the acid/base balance during treatment must be monitored and corrected by changing the pre-dialysis bicarbonate concentration, or by bicarbonate profiling, during treatment. If Citrasate™ is used without reducing anticoagulants, there may be a risk of bleeding and hypercoagulation which in turn may be associated with increased venous access complications.

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