Surefuser™ + Performance data¹

100 mL model ²							
Product reference ³	Nominal volume ⁴	Maximum fill volume ⁵	Flow rate ⁶	Tolerance ⁷	Residual volume ⁸	PCA Bolus volume ⁹	Refill time ¹⁰
P2030/M2345-10P	100 mL	115 mL	2-3-4-5 mL/hr	±15%	2.18 mL	2.0 mL	30 min
P2015/M2345-10P	100 mL	115 mL	2-3-4-5 mL/hr	±15%	2.18 mL	2.0 mL	15 min
300 mL model ²							
Product reference	Nominal volume	Maximum fill volume	Flow rate	Tolerance	Residual volume	PCA Bolus volume	Refill time
P3030/M2345-30P	300 mL	345 mL	2-3-4-5 mL/hr	±15%	4.99 mL	3.0 mL	30 min
P3060/M2345-30P	300 mL	345 mL	2-3-4-5 mL/hr	±15%	4.99 mL	3.0 mL	60 min

- Length of tubing : ¹¹ All types : 900 mm¹²
- 100% filling volume¹³
 - The flow volumes are based on a physiological saline solution with the casing (flow control) at a temperature of 32°C (skin temperature).¹⁴
 - Residual volume was measured when the fluid inside the PCA unit was completely discharged after administration.¹⁵

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SUREFUSER™ + PCA SET

with variable flow regulator
AMBULATORY BALLOON INFUSER

INSTRUCTIONS FOR USE

EN ENGLISH	_____	P1-3
DE GERMAN	DEUTSCH	_____ P4-6
FR FRENCH	FRANÇAIS	_____ P7-9
ES SPANISH	ESPAÑOL	_____ P10-12
IT ITALIAN	ITALIANO	_____ P13-15
NL DUTCH	NEDERLANDS	_____ P16-18
PT PORTUGUESE	PORTUGUÊS	_____ P19-21
EL GREEK	ΕΛΛΗΝΙΚΑ	_____ P22-24
SV SWEDISH	SVENSKA	_____ P25-27
DA DANISH	DANSK	_____ P28-30
NO NORWEGIAN	NORSK	_____ P31-33
FI FINNISH	SUOMI	_____ P34-36
RU RUSSIAN	РУССКИЙ	_____ P37-39
PL POLISH	POLSKI	_____ P40-42
RO ROMANIAN	LIMBA ROMÂNĂ	_____ P43-45
BG BULGARIAN	БЪЛГАРСКИ	_____ P46-48
CS CZECH	ČEŠTINA	_____ P49-51
TR TURKISH	TÜRKÇE	_____ P52-54
HR CROATIAN	HRVATSKI	_____ P55-57
SR SERBIAN	SRPSKI	_____ P58-60
HU HUNGARIAN	MAGYAR	_____ P61-63
LT LITHUANINA	LIETUVIŲ	_____ P64-66
LV LATVIAN	LATVIEŠU	_____ P67-69
SK SLOVAK	SLOVENČINA	_____ P70-72
SL SLOVENIAN	SLOVENŠČINA	_____ P73-75
ET ESTONIAN	EESTI KEEL	_____ P76-78
UK UKRAINIAN	УКРАЇНСЬКА	_____ P79-81

NIPRO

NIPRO CORPORATION
3-9-3, Honjo-Nishi, Kita-ku, Osaka 531-8510, JAPAN

[Manufacturing facility]
NIPRO CORPORATION ODATE FACTORY
8-7, Hanukiyachi, Niida, Odate-shi, Akita, 018-5794 JAPAN

NIPRO MEDICAL EUROPE
Blokhuisstraat 42,
2800 Mechelen, BELGIUM

NIPRO MEDICAL CORPORATION
3150 N.W. 107th Avenue, Miami,
FL 33172, U.S.A.

[Sponsor in Australia]
NIPRO AUSTRALIA PTY LTD
Suite 2.02, Level 2, 657 Pacific Highway,
St Leonards, NSW 2065, AUSTRALIA

[Authorized Representative in Malaysia]
NIPRO MALAYSIA SDN BHD
B-3-2, The Ascent Paradigm, No.1 Jalan
SS7/26A, Kelana Jaya, 47301 Petaling Jaya,
Selangor, MALAYSIA - GC834901272619

Fig. 1 Structure


Fig. 2 Priming Surefuser™ + with Medication


Fig.3 Attachment of Switch Section on Casing and Switch Handle


Fig. 4 Infusion Line Filter


Fig. 5 Priming PCA Unit with Medication


Fig. 6 PCA Unit



Fig. 7 Operation of PCA Unit for Bolus Infusion

* Please note that the shape of parts may change.

Caution

- 1) This is a disposable, single-use-only product that must not be re-sterilized and reused. Discard the product immediately after use.
*Reuse or reprocessing of a single use device may lead to contamination and compromised device function or structural integrity.
- 2) Never fill the balloon reservoir beyond the specified liquid amount, doing so could cause the balloon reservoir to rupture.
- 3) Hold the Surefuser™ + upright with both hands and fill slowly.
- 4) Do not pull excessively on the infusion line; it might cause leakage or the infusion line to detach.
- 5) Air between the balloon reservoir and the filter is expelled by the filter.
- 6) If there is air in the infusion line, the liquid may not flow when the closable clamp is released.
- 7) Make sure there are no bends or twists in the infusion line or the connected catheter. Bends and twists can cause fluctuations in the infusion rate.
- 8) After the balloon reservoir is filled with liquid, use immediately.
- 9) The infusion rate will vary due to changes in the viscosity and density of the medication, temperature and arterial pressure. Please keep these factors in mind when using Surefuser™ +.

SUREFUSER™ + PCA SET

with variable flow regulator

100mL / 300mL

AMBULATORY BALLOON INFUSER

Instructions for Use

Please read these instructions for use carefully before use.

Indications

Surefuser™ + is to be used for continuous drug infusion therapy: post operative pain control, carcinomatous pain control, chemotherapy for cancer etc. All variants have the patient-controlled analgesia (PCA) infuser. Follow instructions by medical professionals.

Contraindications

1. Oil-based medications, i.e., iodized ethyl ester of fatty acid obtained from poppy seed oil, etc., and etoposide medications and fatty emulsion medications should not be used.
2. This product is intended for single use only.

Operating Principle

A rubber balloon reservoir is filled and expanded with medication. The contractile force of the balloon forces out the medication. By running the medication through the casing (flow control) containing a very narrow tube, a specific flow rate can be maintained. The PCA unit contains a reservoir. When used in combination with the PCA unit, the medication that flows out of the balloon reservoir enters the reservoir in the PCA unit. When the bolus infusion button is pressed, the medication will flow out, and the reservoir will be emptied. After the medication completely flows out, the reservoir will be filled with the medication again.

PCA Specifications

- Refill time: approx. 60 min, 30 min, 15min
- Bolus volume: 3.0 mL, 2.0 mL

Method

1) Medication Procedures

1. Refer to the drug insert from the pharmaceutical manufacturer for information about medication procedures.
2. Before beginning patient infusion, make sure you thoroughly understand the infusion characteristics of Surefuser™ + being used.
3. We recommend using Luer lock type syringe.

2) Priming Surefuser™ + with Medication (See Fig. 1-4 below)

1. All operations must be performed under aseptic conditions.
2. Fill the syringe with medication. Completely bleed all air bubbles from the syringe.
3. If a syringe with needle is being used, remove the needle from the syringe.
4. The casing (flow control) is set to the maximum infusion speed when the product is shipped from the factory. Verify the position of the black triangle mark (▲-mark) on the switch. If not set properly, use the switch handle (attached to the back side of the casing) to align the ▲-mark to the maximum infusion speed (See Fig.3). Afterwards, attach the switch handle again to the casing (flow control).
5. Verify correct closing of the closable clamp.
6. Remove the port cap retaining it for later use and fill the balloon reservoir with the medication via the port. Filling operation must be done following the instruction: Place the syringe plunger on a solid surface and place both hands on the syringe for a secure grip. By pressing the syringe barrel slowly down, the solution flows into the balloon reservoir. Pressure must be applied only on the syringe, not on Surefuser™ +. Check the port section for leaks and make sure that there is no damage in the balloon reservoir. The use of a filter is recommended during filling operation. Make sure that the syringe does not separate from the port during filling. (See Fig.2)
7. When the balloon reservoir is filled with the correct liquid volume, disconnect the syringe and close the port with the retained cap.

8. Hold the infusion line filter so that the patient side is at the top. Open the closable clamp on the infusion line to prime the infusion line. After opening the closable clamp, the liquid automatically flows through the product. It is not necessary to remove all of the air from the back side of the filter; air in the back side cannot flow to the patient side.
9. If the liquid does not flow through the product, tap the casing (flow control) with your finger to expel the air.
10. Priming is complete when all air bubbles have been bled from the infusion line and the liquid starts flowing from the connector.
11. When priming is complete, close the infusion line with the closable clamp and replace the cap with the occlusive cap.
12. If necessary, fill out the patient label with the required information and attach the label to the protector.
13. Infusion should be started as soon as possible after the balloon body has been primed.

3) Priming PCA Unit with Medication

1. Close the tube using the closable clamp as shown in the figure below, and then remove the PCA priming port cap. Prepare a 10 mL syringe filled with the medication or its dilute, securely connect it to the female connector of the PCA priming port, and then infuse the medication slowly into the line by the syringe. (See Fig. 5)
2. Fill the infusion line with medication up to the end so that no air remains in it.
3. Remove the syringe from the PCA priming port, and then securely close the PCA priming port cap.
4. Take off the tab on the button area. (See Fig. 6)
NOTE: Using the product with the tab still attached will cause excessive administration of medication inside the balloon. If the tab is removed before use, the product may not have been sterilized well. In such cases, suspend usage of the product immediately.

4) Administering the Medication

1. Use the switch handle to align the ▲-mark to the desired infusion speed. Make sure this is properly done. The flow of medication will stop if the ▲-mark is at an intermediate position.
2. Verify that there are no air bubbles in the infusion line and then attach the connector to the patient's line.
3. Use adhesive tape, etc., to attach the casing (flow control) securely to the patient's skin.
NOTE: If the casing (flow control) is not attached securely to the patient's skin, the medication flow rate may vary from the one desired.

4. Open the closable clamp and start the medication infusion.
5. The scale on the balloon body is only an indication of the volume in the balloon.

5) Changing the infusion Speed (See Fig.3)

1. Close the closable clamp.
2. Remove the casing (flow control) attached to the patient's skin.
3. Use the switch handle to change infusion speed to the desired one.
4. Continue with step 3 from 4) Administering the Medication.

6) Operation of PCA Unit for Bolus Infusion

1. Press the button with constant force until it is fully depressed. (See Fig. 7)
NOTE: Pressing the button but releasing it before it is fully depressed will reduce the infusion amount.
2. Confirm that nothing prevents the return of the button to its original state for the duration of the refill time.
NOTE: Operating the PCA unit before the refill time elapses will reduce the infusion amount.
3. For the next PCA operation, follow the procedure in Operation of the PCA Unit for Bolus Infusion.

- 10) The stated infusion rate for Surefuser™ + is based on using a physiological saline solution with the casing (flow control) at a temperature of 32°C (skin temperature). In order to maintain consistent medication viscosity, the casing (flow control) should be attached securely to the patient's skin using adhesive tape, etc. Infusion rate accuracy is ±15%.
- 11) The specified priming time for Surefuser™ + is based on the use of physiological saline solution at normal temperatures. Additional priming time will be needed when using high-viscosity medications or when using Surefuser™ + at low temperatures.
- 12) Priming should be done at the maximum infusion speed. Priming at anything other than the maximum infusion speed may result in incomplete priming.
- 13) Any modifications to the infusion speed must be made by a physician or qualified personnel. The switch handle should only be entrusted to such persons. Make sure the ▲-mark is properly aligned to the desired infusion speed.
- 14) Before closing the closable clamp, verify that the infusion line is in the center of the closable clamp. If the closable clamp does not clip the infusion line properly, the medication flow will not be stopped.
- 15) Oil-based medications, i.e., iodized ethyl ester of fatty acid obtained from poppy seed oil, etc., and etoposide medications can cause the balloon reservoir to rupture and should not be used with Surefuser™ +.
- 16) Do not use alcohol-containing medication on the filter for disinfection. It could damage its hydrophobic character and could lead to leakage.
- 17) During usage, check the condition of the Surefuser™ + and infusion line regularly for blood clots and leakage. Pay particular attention to each connection part for breakage, leakage or weak connections.
- 18) Using fatty emulsion medications can cause the infusion line filter to become clogged, so the use of such medications should be avoided.
- 19) Refer to the pharmaceutical manufacturer's drug insert for instructions on medication use and quantity.
- 20) If any abnormalities, such as balloon reservoir rupture, medication leaking into the protector, etc., are noticed during use, discontinue use immediately.
- 21) Never use a product if its package is damaged or if any product damage or Irregularities are noticed.
- 22) Some products/medication used can cause the connecting part to be loosened or to be broken. (e.g. the use of alcohol on the connecting parts can cause cracks.)

- 23) Store the product in a cool, dry place away from direct sunlight and excess humidity.
- 24) Select the administration volume and the flow rate of the medication in use as specified for that medication.
- 25) During filling the PCA unit, the button gradually returns to its original position, and therefore ensure not to obstruct the button when this happens. If the return of the button to its original state is obstructed, the medication may not fill the PCA unit.
- 26) Do not press the button before the refill time is finished. The required PCA bolus volume cannot be acquired.
- 27) After operating the PCA unit, do not operate it again until the medication has filled the PCA unit.
- 28) Dispose of the product in an approved biohazard container as per facility protocol.

Guarantee

- 1) Surefuser™ + is manufactured under strict quality control and quality is assured. We will not be responsible, however, for the injury to a patient or any person or for the damage to any object that is attributed to transport, storage, and operation in your institution.
- 2) If a patient or any person is injured or any object is damaged by the use of Surefuser™ +, we will not be responsible for the injury or damage unless we are clearly identified as being at fault.
- 3) If a patient or any person is injured or any object is damaged by the reuse of Surefuser™ +, we will not be responsible for the injury or damage of any nature.
- 4) We will not be responsible for any injury or damage caused by use of Surefuser™ + after the expiry date mentioned on the packages.

Symbols used for labeling

	Do not re-use
	See Instructions for Use
	Sterilized using ethylene oxide
	Expiry date
	Lot number
	Manufacturer
	Reference number
	Authorized representative in the European Community
	Fragile, handle with care
	Do not use if package is damaged
	Non pyrogenic

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