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en: Instructions for use Bicarbonate Haemodialysis Concentrate (B-Concentrate)

Indications: Preparation of dialysis fluid for the treatment of patients with acute or chronic renal failure when haemodialysis or related therapies are prescribed by a physician. **Intended use:** Liquid bicarbonate concentrate for bicarbonate haemodialysis. Use only in combination with an acid haemodialysis concentrate (A-concentrate) in the given dilution (see relevant machine operator's manual). **Side effects:** Bicarbonate concentrate itself does not have known side effects. For side effects of dialysis treatment see acid haemodialysis concentrate. **Contraindications:** Bicarbonate concentrate itself does not have known contraindications. Patients who are hypersensitive to bicarbonate or any components used in the bicarbonate concentrate should not be dialyzed with this product. **Warnings:** Overmixing of bicarbonate concentrate must be avoided. The concentration of the ready-to-use dialysis fluid must be precisely monitored. Before dialysis treatment and each time new concentrate is used with the machine check the conductivity and pH-level of the dialysis

MTN Neubrandenburg GmbH
Gustav-Kirchhoff-Str. 2
17033 Neubrandenburg, Germany
www.nipro-group.com

CE 2797
UK 0086
CA

UK Responsible Person:
Nipro Medical UK Ltd.
Units 12-14 South Point,
Ensign Way, Hamble,
Hampshire
SO31 4RF, United Kingdom

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fluid. Bicarbonate and acid concentrate should only be diluted with water of suitable quality according to ISO 23500-3 immediately before use. Opening the containers can initiate growth of microorganisms, consume within 24 hours. Do not use after the expiration date. Only use undamaged containers and clear colorless liquid. **Note:** The concentration of the ready-to-use dialysis fluid is stated on the label of the corresponding acid haemodialysis concentrate (A-concentrate). Follow the dialysis machine manufacturer's instruction. Depending on the volume, the concentrate can be used for more than one dialysis treatment. Keep container closed when not in use. Concentrate is endotoxin free. For safe application and disposal, see manufacturer's safety data sheet. Any serious incident that has occurred in connection with the product must be reported to the manufacturer and the competent authority of the state in which the user is established.

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