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
EN: ACID CONCENTRATE FOR HEMODIALYSIS

Important:

Store at a temperature no lower than 4°C. Pyrogen-free product. Single use only.

Intended purpose:

Preparation of dialysis fluid used in hemodialysis and related therapies.

Indication:

Treatment of patients with acute or chronic renal failure when hemodialysis is prescribed by a physician.

Clinical Benefits:

Facilitate the clearance of blood from uremic toxins that are able to pass through the dialysis membrane and balance out the electrolyte disbalance.

Contraindication:

There are no special contraindications for the use of Niprobag. Patients who are hypersensitive to any components used in the Niprobag should not be dialyzed with this product.

Side-effects:

There are no side-effects reported for Niprobag.

Instructions:

Acid concentrate with dilution 1:44 to be used in combination with powder bicarbonate concentrate and hemodialysis machines. Niprobag is compatible with Nipro hemodialysis machines. Before using, consult the hemodialysis machine manual and follow the instructions for connecting and disconnecting.

Remove the cap. Connect the bag tube to the acid connector on the machine following the instructions of the hemodialysis machine manufacturer. Break the bag connector ("luer lock") that is inserted inside the tube.

Dilution is based on the following proportion: See the product label. B is the bicarbonate concentrate and H2O is deionized water, complying with EP, ISO 23500-5, and applicable to local quality requirements.

Warnings:

To be used by healthcare professionals. Home dialysis is excluded. To be used under the physician's order. Make sure the composition corresponds with the physician's prescription. Do not use if package is damaged. The concentrated solution must be diluted immediately before use. Any quantity of unused solution must be discarded as neither its properties nor its final mixture composition can be guaranteed, and therefore the best treatment for the patient cannot be assured. See table below for an indication of expected dialysis fluid volume obtained from each package volume. Do not use after the expiry date stated on the package. Any serious incident that has occurred in relation to the product should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is based. Dispose of the product according to local and national regulations.

SIZE	DILUTION FACTOR	DIALYSIS FLUID VOLUME
3.5L	1:44	157.5L
3.8L	1:44	171.0L
4.7L	1:44	211.5L

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