GOODTEC Y-CONNECTOR SET

EN	English	English	2
DA	Danish	Dansk	10
DE	German	Deutsch	13
ES	Spanish	Español	16
	Spanish (LATAM)		19
FI	Finnish	Suomi	22
FR	French	Français	25
EL	Greek	Ελληνικά	28
IT	Italian	Italiano	31
NL	Dutch	Nederlands	34
NO	Norwegian	Norsk	37
PT	Portuguese	Português	40
	Portuguese (LATAM)		43
SV	Swedish	Svenska	46
LT	Lithuanian	Lietuvių K.	49
RU	Russian	Русский	52
CS	Czech	Česky	55
HU	Hungarian	Magyar	58
PL	Polish	Polski	61
SK	Slovak	Slovensky	64
ET	Estonian	Eesti Keel	67
LV	Latvian	Latviešu	70
TR	Turkish	Türkçe	73
	Labeling Symbol Definitions		76

Manufacturer

Goodman Co., Ltd.

5F KDX Nagoya Sakae Building,
4-5-3 Sakae, Nakaku, Nagoya, Aichi
460-0008 Japan

Manufacturing site

Goodman Co., Ltd.

Goodman Research Center

276-1 Idogane-cho, Seto, Aichi
489-0976 Japan

Authorized representative in EU

Goodman Medical Ireland Ltd

Mervue Business Park
Galway H91 H9CK Ireland

UK Responsible person

Emergo Consulting (UK) Limited

c/o Cr360 – UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ United Kingdom

CH authorised representative

MedEnvoy Switzerland

Gothardstrasse 28
6302 Zug Switzerland

21 May 2024/Rev 2/GM25/L01750003X

ENGLISH

Warning

1. This product is sterilized with ethylene oxide gas (EOG), is not reusable and must not be re-sterilized. Re-sterilization could lead to infection of patient or degradation of the material characteristics of the device.
2. Only use this device after air has been completely removed from lumen. Complications such as air embolism may occur.
3. The fixed valve should be in the closed position when negative pressure is applied to the Y-connector to prevent air from being drawn into the device. Complications such as air embolism may occur.
4. Both hemostatic valve and fixed valve should be completely open when inserting or removing catheters or other medical devices. Damage may occur to Y-connector and/or catheter, etc.
5. Ensure that air has been entirely removed from catheter, etc. prior to insertion into Y-connector. Complications such as air embolism may occur.
6. If any resistance is felt when opening the fixed valve, immediately stop rotating valve. Damage to Y-connector may occur.
7. Not to be used with products comprising of organic solvent, fat emulsion, or oil-based components. It may incur damage to product.
8. Injection pressures of greater than 3448kPa (500psi) should not be performed when infusing contrast media. Excessive pressure may cause damage to product.
9. Do not connect the male connector of the included extension tube to anything other than the Y connector side port.

Intended patient population

Patients who undergo treatment or examination by catheterizations.

Clinical benefit

The device assists users in performing an optimal series of complicated medical operations performed by hand while reducing patient burden by leak-proofing, and contributes to management of therapeutic treatment through side-port monitoring.

Shape/Construction

Details are shown in the table below.

List of Models

Set model	Type of Y-connector	Included accessories		
		Inserter	Torque device	Extension tube
YOK0A	OKAY II	●	●	
YOK0B	OKAY II			
YOK0E	OKAY II	●	●	●
YOKS0A	OKAY II Short	●	●	
YOKS0B	OKAY II Short			
YOKS0E	OKAY II Short	●	●	●
YSS0A	OKAY II Super Short	●	●	
YSS0B	OKAY II Super Short			
YSS0E	OKAY II Super Short	●	●	●

1. Product Summary

(1) Y-connector

The Y-connector is attached to a device such as a guiding catheter for the purposes of reducing blood leakage and for fixation of devices inserted into the vasculature. Simultaneously, a channel is provided for obtaining blood pressure measurement or the introduction of drug solution/contrast media.

There are three types of Y-connectors with different lengths.

(2) Inserter

The inserter facilitates smooth introduction of guide wires into the Y-connector without causing damage to the tip shape.

(3) Torque device

The torque device is attached appropriately to the proximal section of the guide wire in order to assist in delivery of guide wire to target lesion. The torque device provides for easy gripping when rotating or manipulating the guide wire.


(4) Extension tube

The extension tube is attached to the sideport and other devices used for passing contrast media and other solutions.


2. Construction

(1) GOODTEC Y-connector


A. OKAY II

B. OKAY II Short Type

C. OKAY II Super Short Type


(2) Inserter

(3) Torque device

(4) Extension tube

3. Principles of operation

(1) Hemostatic valve

Pushing the opener in the direction of the hemostatic valve opens the device and pulling it in the other direction closes it. Pushing the opener and turning in a clockwise direction locks the device in an open position. While the hemostatic valve is closed, blood loss during catheter manipulation is kept to a minimum.

(2) Fixed valve

Rotating the screw in a right direction closes the valve while rotating in a left direction opens it. With the fixed valve in a closed position, catheter

positioning can be fixed and introduction of contrast media or drug solution is possible.

Purpose

This product is designed to work in conjunction with guiding catheters, etc. to reduce blood loss while assisting in manipulation of catheters, etc., injection of contrast solution from a side port, injection of medicine or saline solution, and measurement of blood pressure.

Specifications

1. Inserter tensile strength: No detachment when a tensile force of 35N at 10N/sec is applied for a 15 second duration
2. Torque device guide wire fixation strength: >9.8N
3. Main Branch ID: 3.33mm/10Fr (The ID of the valve segments are less than 3.33mm/10Fr.)
4. Female connectors conform to ISO80369-7
Rotator and extension tube male connectors conform to ISO 594-1,-2
The clinical risk is determined to fall within the acceptable parameters.

Method of use

1. All equipment and devices to be used should be carefully inspected and proper functionality confirmed prior to use.
2. A line is attached to the Y-connector side port for the purpose of obtaining blood pressure measurement, introduction of contrast media, or medication.
3. In order to eliminate air from within the Y-connector, the hemostatic valve and fixed valve should be in the open position, a finger placed across the rotator opening and heparinized saline flushed through the side port.
4. In order to avoid the introduction of air after emission, flushing should be performed while closing the hemostatic valve, and the Y-connector should be filled with heparinized saline.
5. The Y-connector rotator section is attached to the guiding catheter hub.
6. When introducing a guide wire into the Y-connector, the tip of the guide wire should be placed within the inserter, the hemostatic valve opened and the inserter placed deep within the Y-connector.
7. After insertion, the inserter should be removed while leaving the guide wire in place and the hemostatic valve closed.

8. The torque device is attached to an appropriate position on the guide wire, the guide wire advanced to the target lesion and the torque device removed from the guide wire.
9. For introducing catheters after being placed on the guide wire, the Y-connector hemostatic valve should be opened and the catheter inserted. After insertion of the main catheter component, the hemostatic valve should be closed and catheter delivered to target location.
10. When catheter, etc. reaches target location, the fixed valve may be closed to fixate the catheter.

Precautions for use

1. Precautions prior to use

- (1) Be sure to reference IFUs accompanying all medical devices or products to be used jointly with this device.
- (2) Confirm the compatibility of this device and other medical devices to be used in regards to the procedure.
- (3) If the package or product is damaged or contaminated, do not use the product.
- (4) All device use should be undertaken in a sterile environment.
- (5) One of the components contains the following substance defined as Carcinogenicity; Category 1B and Reproductive toxicity; Category 1B at a concentration exceeding 0.1% by weight.

Cobalt: CAS No. 7440-48-4

Current scientific evidence supports that exposure to cobalt-containing alloys in medical devices does not present a risk of carcinogenicity or reproductive toxicity.

2. Precautions during use

- (1) The hemostatic valve is designed to facilitate catheter delivery, and will not completely eliminate blood loss, but is instead designed to minimize it. Conditions of use will determine volume of blood loss.
- (2) Introduction of contrast media or drug solution through the Y-connector should be performed with the fixed valve in a closed position and within the provisions of infusion pressure. Exceeding infusion pressure may result in the leakage of contrast media or drug solution.

- (3) Do not over tighten the fixed valve when fixing catheters into position. Damage to catheter may occur.
- (4) If resistance is felt when introducing catheters through the Y-connector, immediately cease insertion and confirm the cause of the resistance. Damage to catheter may occur.
- (5) The hemostatic valve of the OKAY II type product is not constructed to close automatically. The opener must be pulled to the proper position by hand when closing the hemostatic valve.
- (6) The inner lumen and outer casing of the Y-connector should be regularly flushed with heparinized saline. Contrast media and/or blood may adhere to device, affecting closing functionality of the hemostatic valve and/or the fixed valve.
- (7) Should the hemostatic valve closing function be affected by contrast media or blood adherence, blood loss can be controlled through manipulation of the fixed valve when inserting catheters.
- (8) When an extension tube is utilized, do not over tighten connection. Damage to y-connector extension tube may result.

3. Precautions post use

Dispose of this product as medical waste, while considering the prevention of infection.

4. Adverse event

- Hypotension/hypertension
- Hemorrhagic complications
- Infection
- Air embolism

Storage method, shelf life and Other

1. Storage method

- (1) Store the product in a room temperature location not exposed to high temperature and humidity or direct sunlight and take proper precautions to ensure product does not contact water.
- (2) Avoid inclinations, vibrations and impacts (including during transportation) and store in a safe, stable environment.

- (3) Do not store near chemicals or in areas where the device may be exposed to gases.

2. Shelf-life

Use this product before the “expiration date” shown on the package label.

PRODUCT WARRANTY DISCLAIMERS AND LIMITATION OF REMEDIES

With respect to any and all Goodman product(s) depicted or described in this publication, there are no express or implied warranties of any kind whatsoever, including, without limitation, any implied warranty of merchantability or fitness for a particular purpose. Goodman Co., Ltd and its subsidiaries (“Goodman”) shall not be liable for any direct, incidental, consequential, or other damages in connection with the use, re-use, or any other aspect of the relevant product(s) except as expressly provided under applicable law. No person has the authority to bind Goodman to any representation or warranty. Any product or service description or specifications contained in any Goodman printed matter, including this publication, are meant solely to generally describe the relevant product at the time of manufacture and do not constitute any express or other warranties.

Any incidents shall be reported to the manufacturer and the competent authority of your State.

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