



## SYRINGE WITH LUER LOCK, WITHOUT NEEDLE

STERILE, SINGLE-USE SYRINGE FOR INJECTION OR ASPIRATION

Manufactured in Shanghai, China



# SYRINGE WITH LUER LOCK, WITHOUT NEEDLE

STERILE, SINGLE-USE SYRINGE FOR INJECTION OR ASPIRATION

- Transparent barrel for easy recognition of air bubbles
- Smooth plunger motion along entire length of barrel
- Large finger grips for secure handling
- Latex-free, BPA-free, DEHP-free, PVC-free
- Available in centric tip

CE<sub>0123</sub>

Class I Medical Devices

Rule 1

MDD 93/42/EEC

UMDNS 13929

## PRODUCT COMPLIANCE

- CE-marked, Class I Medical Devices, Rule 1, MDD 93/42/EEC, UMDNS 13929
- Complies with the following norms, directives, and regulations:
  - ISO 594-1
  - EN 1707
  - ISO 7886-1
  - ISO 13485
  - ISO 10993-1
  - ISO 10993-4
  - ISO 10993-5
  - ISO 10993-7
  - ISO 11607-1
  - ISO 11607-2
  - EN ISO 13485
  - ISO 15223-1
- Labels contain 15 languages

## MANUFACTURING DETAILS

**Legal manufacturer:** Shanghai Kindly Enterprise Development Group Co., Ltd.

**Country of origin:** China

## STERILIZATION AND SHELF LIFE

### Sterilization

EtO (Ethylene oxide)



### Single-use only




### Shelf life

5 years



## MATERIALS USED



DEHP-free  
BPA-free  
Latex-free  
PVC-free

## PRODUCT RANGE OVERVIEW

Reference code	Description	Size (ml)	Tip type
SY3-1LC-EC	3 part syringe 1 ml luer lock	1	CENTRIC
SY3-2LC-EC	3 part syringe 2 ml luer lock	2	CENTRIC

## PACKAGING DETAILS

Reference	Size	Inner box				Outer box				Net weight	Gross weight
		Quantity	L (mm)	W (mm)	H (mm)	Quantity	L (cm)	W (cm)	H (cm)		
SY3-1LC-EC	1	100	210	165	92	3600	58	44	53	7,9	10,9
SY3-2LC-EC	2	100	420	125	70	3200	58	44	53	12,0	18,0

### Languages

15 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), Finnish (FI), Russian (RU), Turkish (TR), and Polish (PL)





### Transport conditions

Closed and dry

### Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
REF	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
	—	—	✓
	✓	✓	✓

