



## SUREFUSER™ + BAG

EASY TO WEAR BAG FOR SUREFUSER™ +  
ELASTOMERIC INFUSION PUMP



# SUREFUSER™ + BAG

EASY TO WEAR BAG FOR SUREFUSER™+ ELASTOMERIC INFUSION PUMP

- Waterproof
- Thermo insulated
- UV protected
- Lightweight
- Easy to clean
- Latex free
- Available in four sizes
- Makes pump portable to improve mobility
- Easy to wear around waist or shoulder
- Adjustable strap
- Zipper opening for convenient access
- Long version available when Surefuser is equipped with a CSTD

## MATERIALS USED

**Material:** Elastic sponge with UV foil

**Material:** Fabric - Oxford



## PRODUCT RANGE OVERVIEW

Reference code	Size	Intended for
SFS-BAG-L-PL	Large	pumps ranging from 50 - 100 ml
SFS-BAG-XL-PL	X-Large	pumps ranging from 150 - 300 ml
SFS-BAG-L-LONG-PL*	Large Long	pumps ranging from 50 - 100 ml
SFS-BAG-XL-LONG-PL*	X-Large Long	pumps ranging from 150 - 300 ml

\*The longer version of our Surefuser bags is specially designed to fit both the pump & connector when the Surefuser is equipped with a Closed System Transfer Device (CSTD)

## PACKAGING DETAILS

Reference code	Size	Length	Width	Height	Weight	Quantity
Individual	Large	19 cm	5,5 cm	-	50 G	-
	X-Large	19 cm	9 cm	-	60 G	-
	Large Long	23 cm	5,5 cm	-	60 G	-
	X-Large Long	23 cm	9 cm	-	70 G	-
Outer box	Large	60	54	45	21,5 kg	400 pces
	X-Large	60	54	45	16,5 kg	250 pces
	Large Long	60	54	45	16,5 kg	250 pces
	X-Large Long	60	54	45	19 kg	250 pces

Country of origin: Poland

## TRANSPORT AND STORAGE CONDITIONS

Closed and dry

SureFuser+ elastomeric pump not included. For details about SureFuser+, please refer to the corresponding Product Data Sheet.

