



## GM-30N™

INFLATION DEVICE




## Features:

Complying with the most stringent of test parameters, the GM-30N is representative of the emphasis Nipro places on producing reliable, high quality products, while the ergonomic design ensures ease-of-use.

### • Pressure gauge

Maximum pressure reading:  
3.0 MPa (30 atm / 441 psi).

Fluorescent display allows for viewing in darkened conditions.



### • Release lever

Ergonomic design that fits comfortably in the hand to allow quick release with minimum effort.

### • Handle

Designed for easy rotation and handle extraction.



# Product information

Maximum Pressure Reading	3.0 MPa (30 atm / 441 psi)
Maximum Volume Dispensed	25 ml
Accuracy	±3 %FS

# Ordering information

Model No.
GM-30N: Three-way stopcock, turns 180°
GM-30NF: Three(four)-way stopcock, turns 360°



Nipro Corporation Japan is a leading global healthcare company established in 1954. With over 35,000 employees worldwide, Nipro serves the Medical Device, Pharmaceutical, and Pharmaceutical Packaging industries.

Nipro is dedicated to providing interventional cardiologists and vascular surgeons with tools to perform delicate lifesaving and life-enhancing procedures to treat patients with vascular and coronary diseases.

We develop, manufacture, and market cutting-edge technology as well as the newest techniques, innovative solutions, and complex advancements to reshape the vascular space globally. Working hand in hand with physicians, we are focused on providing the best possible care for patients.

Nipro is dedicated to advancing the diagnosis and management of coronary artery disease worldwide as part of our promise to help people.

**LIVE LONGER. LIVE BETTER.**

