



DOSESAVER™ ND

1 ML SYRINGE WITH FIXED NEEDLE

DEAD VOLUME OF ONLY 0,002 ML



Dosesaver™ ND

1 ml syringe with fixed needle

DEAD VOLUME OF ONLY 0,002 ML

- Syringe tip space eliminated by fixed needle
- Available in 1 ml volume and 25 G 25 mm needle
- Three-beveled lancet cut needle to minimize discomfort
- Grey gasket instead of the regular black gasket for better visibility of the scale
- Reduces risk of over-sampling
- Reduces drug waste due to low dead volume

CE₀₁₂₃
Class IIa
Rule 6
Annex IX

Product Compliance

- Complies with the following norms, directives, and regulations:
 - ISO 6009:2016
 - EN ISO 10993-1:2020
 - EN ISO 11135:2014/A1:2019
 - ISO 7864:2016
 - EN ISO 10993-4:2017
 - EN ISO 11607-1:2020
 - EN ISO 15223-1:2016
 - EN ISO 10993-10:2013
 - EN ISO 11607-2:2020
 - EN ISO 80369-7:2021
 - EN ISO 10993-11:2018
 - ISO 80369-7:2021

Manufacturing Details

Legal manufacturer: Nipro Corporation

Country of origin: Thailand

Sterilization and Shelf Life

Sterilization

EtO (Ethylene oxide)



Single-use only




Shelf life

5 years



Materials Used



PRODUCT RANGE OVERVIEW

| Reference code | Description | Dead space volume (ml) |
|------------------|---|------------------------|
| SY3-ND-1-2525-ET | Dosesaver ND Syringe ND 1 ml fixed needle 25 G 25mm | 0,002 ml |

Packaging Details

| Individual unit | | | | Inner box | | | | Outer box | | | | |
|-----------------|------------|-------------|------------|-------------|------------|-------------|----------|-------------|------------|-------------|-------------------|----------|
| Length (mm) | Width (mm) | Height (mm) | Weight (g) | Length (mm) | Width (mm) | Height (mm) | Quantity | Length (mm) | Width (mm) | Height (mm) | Gross weight (kg) | Quantity |
| 166 | 45 | 9 | 5,49 | 232 | 497 | 372 | 100 | 580 | 380 | 250 | 5,49 | 1000 |

Languages

26 languages on inner box:

Arabic (AR), Bulgarian (BG), Czech (CS), Danish (DA), English (EN), German (DE), Greek (EL), Spanish (ES), Estonian (ET), Finnish (FI), French (FR), Croatian (HR), Hungarian (HU), Italian (IT), Lithuanian (LT), Latvian (LV), Dutch (NL), Norwegian (NO), Polish (PL), Portuguese (PT), Romanian (RO), Russian (RU), Slovakian (SK), Slovenian (SL), Serbian (SR), Swedish (SV), Turkish (TR)

Storage conditions



Do not store at extreme temperature and humidity.






Not made with
latex, DEHP, PVC or BPA

Transport conditions

Closed and dry.

Label details

| | Individual unit | Inner box | Outer box |
|---|-----------------|-----------|-----------|
| Product name | ✓ | ✓ | ✓ |
| Product size | ✓ | ✓ | ✓ |
| Quantity (pieces) | — | ✓ | ✓ |
| REF | — | ✓ | ✓ |
| LOT | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
| STERILE EO | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |

| | Individual unit | Inner box | Outer box |
|---|-----------------|-----------|-----------|
| CE | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
| EC REP | ✓ | ✓ | ✓ |
|  | — | ✓ | ✓ |
| Made in Thailand | ✓ | ✓ | ✓ |

