



DOSESAVER™ LD

1 ML SYRINGE WITHOUT NEEDLE

LOW DEAD VOLUME OF ABOUT 0,062 ML



Dosesaver™ LD

1 ml syringe without needle

LOW DEAD VOLUME OF ABOUT 0,062 ML

- Available in 1 ml volume
- Grey gasket instead of the regular black gasket for better visibility of the scale
- Reduces risk of over-sampling
- Reduces drug waste

CE₀₁₂₃

Class I

Medical Devices

Rule 1

MDD 93/42/EEC

UMDNS 13929

Product Compliance

- Complies with the following norms, directives, and regulations:
 - ISO 7864:2016
 - EN ISO 80369-7:2021
 - EN ISO 11607-1:2019
 - EN ISO 11607-2:2019
 - ISO 80369-7:2021
 - ISO 7886-1: 2017
 - ISO 80369-7: 2021
 - ISO 11607-1: 2019
 - ISO 11607-2: 2019
 - ISO 10993-7: 2008/Amd 1: 2019

Manufacturing Details

Legal manufacturer: Nipro Corporation

Country of origin: Thailand

Sterilization and Shelf Life

Sterilization

EtO (Ethylene oxide)



Single-use only



Shelf life

5 years



Materials Used



PRODUCT RANGE OVERVIEW

Reference code	Description
SY3-LD-1SC-ET	LD 3 part syr 1 ml W/O needle slip

Not made with
latex, DEHP, PVC or BPA

Packaging Details

Individual unit			Inner box				Outer box				
Length (mm)	Width (mm)	Weight (g)	Length (mm)	Width (mm)	Height (mm)	Quantity	Length (mm)	Width (mm)	Height (mm)	Weight (kg)	Quantity
140	96	4,62	202	135	91	100	229	482	337	4,62	1000

Languages

27 languages on all levels of packaging:

Arabic (AR), Bulgarian (BG), Czech (CS), Danish (DA), English (EN), German (DE), Greek (EL), Spanish (ES), Estonian (ET), Finnish (FI), French (FR), Croatian (HR), Hungarian (HU), Italian (IT), Lithuanian (LT), Latvian (LV), Dutch (NL), Norwegian (NO), Polish (PL), Portuguese (PT), Romanian (RO), Russian (RU), Slovakian (SK), Slovenian (SL), Serbian (SR), Swedish (SV), Turkish (TR)

Storage conditions

Do not store at extreme temperature and humidity.

Transport conditions

Closed and dry.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
REF	—	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
	✓	✓	✓

	Individual unit	Inner box	Outer box
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
Made in Thailand	✓	✓	✓

