



DOSESAVER™ LD

1 ML SYRINGE WITH DETACHABLE NEEDLE

LOW DEAD VOLUME OF ABOUT 0,062 ML



Dosesaver™ LD

1 ml syringe with detachable needle

FOR THE INJECTION OF MEDICATION WITH A LOW DEAD VOLUME OF ABOUT 0,062 ML

- Available in 1 ml and several needle combinations
- Three-beveled lancet cut needle to minimize discomfort
- Grey gasket instead of the regular black gasket for better visibility of the scale
- Reduces risk of over-sampling
- Reduces drug waste

CE 0123

Class IIa Medical Devices

Rule 6

MDD 93/42/EEC

UMDNS 13940

Product Compliance

- Complies with the following norms, directives, and regulations:
 - ISO 6009:2016
 - EN ISO 11135:2014/A1:2019
 - ISO 80369-7:2021
 - ISO 7864:2016
 - EN ISO 11607-1:2019
 - EN ISO 80369-7:2021
 - EN ISO 11607-2:2019

Manufacturing Details

Legal manufacturer: Nipro Corporation

Country of origin: Thailand

Sterilization and Shelf Life

Sterilization

EtO (Ethylene oxide)

STERILE EO

Single-use only




Shelf life

5 years



Materials Used



PRODUCT RANGE OVERVIEW

Reference code	Description
SY3-LD-1SC-2325-ET	LD 3 part syr 1 ml with needle 23 G 25 mm
SY3-LD-1SC-2525-ET	LD 3 part syr 1 ml with needle 25 G 25 mm
SY3-LD-1SC-2538-ET	LD 3 part syr 1 ml with needle 25 G 38 mm
SY3-LD-1SC-2516-ET	LD 3 part syr 1 ml with needle 25 G 16 mm
SY3-LD-1SC-2613-ET	LD 3 part syr 1 ml with needle 26 G 13 mm
SY3-LD-1SC-2713-ET	LD 3 part syr 1 ml with needle 27 G 13 mm

Not made with
latex, DEHP, PVC or BPA

Packaging Details

Individual unit				Inner box				Outer box				
Length (mm)	Width (mm)	Height (mm)	Weight (g)	Length (mm)	Width (mm)	Height (mm)	Quantity	Length (mm)	Width (mm)	Height (mm)	Weight (kg)	Quantity
180	96	—	5.43	205	175	88	100	497	372	232	5,93	1000

Languages

27 languages on inner box:

Arabic (AR), Bulgarian (BG), Czech (CS), Danish (DA), English (EN), German (DE), Greek (EL), Spanish (ES), Estonian (ET), Finnish (FI), French (FR), Croatian (HR), Hungarian (HU), Italian (IT), Lithuanian (LT), Latvian (LV), Dutch (NL), Norwegian (NO), Polish (PL), Portuguese (PT), Romanian (RO), Russian (RU), Slovakian (SK), Slovenian (SL), Serbian (SR), Swedish (SV), Turkish (TR)

Storage conditions

Do not store at extreme temperature and humidity.

Transport conditions

Closed and dry.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
REF	—	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
	✓	✓	✓

	Individual unit	Inner box	Outer box
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
Made in Thailand	✓	✓	✓

