



2-PART SYRINGE

STERILE, SINGLE-USE FOR IMMEDIATE
INJECTION OR BLOOD COLLECTION



2-PART SYRINGE

STERILE, SINGLE-USE FOR IMMEDIATE INJECTION OR BLOOD COLLECTION

- Clear barrel for precise volume reading
- Smooth plunger motion along entire length of barrel
- Large finger grips for secure handling
- Latex-free, BPA-free, DEHP-free, PVC-free
- Available in centric or eccentric tip
- For use immediately after filling; not intended to contain medicine for an extended period of time

CE₀₁₂₃
Medical Devices
Class I
UMDNS: 13929
Rule 1

PRODUCT COMPLIANCE

- CE-marked, Class I Medical Devices, Rule 1, UMDNS: 13929
- Complies with the following norms, directives, and regulations:
 - EN 556-1
 - EN 1707
 - ISO 7886-1
 - EN ISO 10993-1
 - EN ISO 10993-4
 - EN ISO 10993-5
 - EN ISO 10993-7
 - EN ISO 10993-10
 - EN ISO 10993-11
 - EN ISO 11135
 - EN ISO 11737-2
 - EN ISO 11138-1
 - EN ISO 11138-2
 - EN ISO 11607-1
 - EN ISO 11607-2
 - EN ISO 11737-1
 - EN ISO 13485
 - ISO 14644-1
 - ISO 14698-1
 - EN ISO 14971
 - EN 15986
 - EN 20594-1
 - EN 62366
- Labels contain 26 languages

MANUFACTURING DETAILS

Legal manufacturer: Guangdong Kindly Medical Device Group Co., Ltd.

Country of origin: China

STERILIZATION AND SHELF LIFE

Sterilization

EtO (Ethylene oxide)



Single-use only



Shelf life

5 years



MATERIALS USED



Barrel



Material: Ethylene-Polypropylene, Copolymer

Plunger

Material: Polyethylene (PE)

DEHP-free
BPA-free
Latex-free
PVC-free

SCALING



PRODUCT RANGE OVERVIEW

Reference code	Description	Volume	Type	Scaling	Max. graduation	Dead space
SY2-2SC-EC	2 part syringe 2 ml Luer Slip	2 ml	CENTRIC	0,1 ml	2,4 ml	≤ 0,055 ml
SY2-5SC-EC	2 part syringe 5 ml Luer Slip	5 ml	CENTRIC	0,2 ml	5,8 ml	≤ 0,061 ml
SY2-5ESC-EC	2 part syringe 5 ml Luer Slip	5 ml	ECCENTRIC	0,2 ml	5,8 ml	≤ 0,061 ml
SY2-10ESC-EC	2 part syringe 10 ml Luer Slip	10 ml	ECCENTRIC	0,5 ml	11 ml	≤ 0,061 ml
SY2-20ESC-EC	2 part syringe 20 ml Luer Slip	20 ml	ECCENTRIC	1 ml	24 ml	≤ 0,081 ml

PACKAGING DETAILS

Box	Quantity (pcs)		Outer box				
	Inner box	Outer box	Net weight (kg)	Gross weight (kg)	Length (mm)	Width (mm)	Height (mm)
SY2-2SC-EC	100	2500	6,65	9,90	600	405	375
SY2-5SC-EC	100	2000	8,54	12,00	605	400	465
SY2-5ESC-EC	100	2000	8,54	12,00	605	400	465
SY2-10ESC-EC	100	1200	8,39	11,00	395	405	605
SY2-20ESC-EC	100	800	9,26	11,85	605	410	465

Languages

26 languages on all levels of packaging:

English (EN), Bulgarian (BG), Czech (CS), Danish (DA), German (DE), Greek (EL), Spanish (ES), Estonian (ET), Finnish (FI), French (FR), Croatian (HR), Hungarian (HU), Italian (IT), Lithuanian (LT), Latvian (LV), Dutch (NL), Norwegian (NO), Polish (PL), Portuguese (PT), Romanian (RO), Russian (RU), Slovakian (SK), Slovenian (SL), Serbian (SR), Swedish (SV), Turkish (TR)

Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
REF	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
	—	✓	✓
	✓	✓	✓
	✓	✓	✓

