



## WING CATH

STERILE, SINGLE-USE WINGED CATHETER FOR SEVERAL HOURS TO DAYS INTRAVENOUS INFUSION OR BLOOD SAMPLING



# WING CATH

STERILE, SINGLE-USE WINGED CATHETER FOR SEVERAL HOURS TO DAYS INTRAVENOUS INFUSION OR BLOOD SAMPLING

- Radiopaque
- Injection port with universal luer fitting
- Flexible wings for easy grip
- Color-coded catheter port to indicate gauge
- Gamma-ray sterilized
- Siliconization of needle-tip to minimize discomfort
- Easy blood visualization in the transparent flashback chamber to facilitate insertion
- Not made with latex or PVC
- For use by healthcare professionals only

CE<sub>0123</sub>

Class IIa Medical Devices  
Rule 7  
MDD 93/42/EEC  
UMDNS 10727

## PRODUCT COMPLIANCE

- CE-marked, Class IIa Medical Devices, Rule 7, MDD 93/42/EEC, UMDNS 10727
- Complies with the following norms, directives, and regulations:
  - ISO 594-1
  - EN ISO 10555-1
  - ISO 15223-1
  - ISO 594-2
  - EN ISO 13485
- Labels contain 12 languages

## MANUFACTURING DETAILS

**Legal manufacturer:** Nipro India Corporation PVT. LTD.

**Country of origin:** India

## STERILIZATION AND SHELF LIFE

### Sterilization

Gamma-ray



### Single-use only



### Indicator sticker

Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a yellow sticker that turns red when sterilized.

### Shelf life


5 years



# MATERIALS USED

**Not made with Latex or PVC**

**Lubricant Material: Silicone**



# PRODUCT RANGE OVERVIEW

Reference code	Description	Purpose	MEAN Flow rate	Color
NIC-1650-PWF-GIN2	IV Cath 16 Gx2" (16x50 mm) ETFE	IV	180 (ml/min)	Grey
NIC-1832-PWF-GIN2	IV Cath 18 Gx1 1/4" (18x32 mm) ETFE	IV	90 (ml/min)	Green
NIC-1850-PWF-GIN2	IV Cath 18 Gx2" (18x50 mm) ETFE	IV	90 (ml/min)	Green
NIC-2025-PWF-GIN2	IV Cath 20 Gx1" (20x25 mm) ETFE	IV	61 (ml/min)	Purple
NIC-2032-PWF-GIN2	IV Cath 20 Gx1 1/4" (20x32 mm) ETFE	IV	61 (ml/min)	Purple
NIC-2050-PWF-GIN2	IV Cath 20 Gx2" (20x50 mm) ETFE	IV	61 (ml/min)	Purple
NIC-2225-PWF-GIN2	IV Cath 22 Gx1" (22x25 mm) ETFE	IV	36 (ml/min)	Blue
NIC-2232-PWF-GIN2	IV Cath 22 Gx1 1/4" (22x32 mm) ETFE	IV	36 (ml/min)	Blue
NIC-2419-PWF-GIN2	IV Cath 24 Gx3/4" (24x19 mm) ETFE	IV	23 (ml/min)	Yellow
NIC-2425-PWF-GIN2	IV Cath 24 Gx1" (24x25 mm) ETFE	IV	23 (ml/min)	Yellow

# PACKAGING DETAILS

Reference	Gauge	Length	Inner box					Outer box				
			Quantity	L (mm)	W (mm)	H (mm)	Weight (g)	Quantity	L (cm)	W (cm)	H (cm)	Weight (kg)
NIC-1650-PWF-GIN2	16	50	100	168	136	150	601	1000	698	320	178	6,880
NIC-1832-PWF-GIN2	18	32	100	168	136	150	573	1000	698	320	178	6,660
NIC-1850-PWF-GIN2	18	50	100	168	136	150	573	1000	698	320	178	6,660
NIC-2025-PWF-GIN2	20	25	100	168	136	150	550	1000	698	320	178	6,640
NIC-2032-PWF-GIN2	20	32	100	168	136	150	550	1000	698	320	178	6,640
NIC-2050-PWF-GIN2	20	50	100	168	136	150	550	1000	698	320	178	6,640
NIC-2225-PWF-GIN2	22	25	100	168	136	150	519	1000	698	320	178	6,120
NIC-2232-PWF-GIN2	22	32	100	168	136	150	519	1000	698	320	178	6,120
NIC-2419-PWF-GIN2	24	19	100	168	136	150	515	1000	698	320	178	6,080
NIC-2425-PWF-GIN2	24	25	100	168	136	150	515	1000	698	320	178	6,080

## Languages

12 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Polish (PL), Portuguese (PT), Finnish (FI), Swedish (SV), Danish (DA), Romanian (RO), Bulgarian (BG)














## Transport conditions

Closed and dry

## Storage conditions

Open the packaging only immediately before use to guarantee sterility.

## Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	—	—
	—	—	✓
	✓	—	—

