



## 3-PART SYRINGE WITHOUT NEEDLE

STERILE, SINGLE-USE SYRINGE FOR INJECTION OR ASPIRATION



# 3-PART SYRINGE WITHOUT NEEDLE

STERILE, SINGLE-USE SYRINGE FOR INJECTION OR ASPIRATION

- Transparent barrel for easy recognition of air bubbles
- Smooth plunger motion along the entire length of the barrel
- Large finger grips for secure handling
- Not made with latex or PVC
- Available in centric or eccentric tip
- Available in luer lock or luer slip

CE<sub>0123</sub>

Class I Medical Devices

Rule 2

MDD 93/42/EEC

UMDNS 13929

## PRODUCT COMPLIANCE

- CE-marked, Class I Medical Devices, Rule 2, MDD 93/42/EEC, UMDNS 13929
- Complies with the following norms, directives, and regulations:
  - ISO 594-1
  - ISO 7864
  - ISO 7886-1
  - ISO 9626
  - ISO 11607-1
  - EN ISO 13485
  - EN ISO 15223-1
- Labels contain 15 languages

## MANUFACTURING DETAILS

**Legal manufacturer:** Nipro India Corporation PVT. LTD.

**Country of origin:** India

## STERILIZATION AND SHELF LIFE

### Sterilization

Gamma-ray



### Single-use only



### Indicator sticker


Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a yellow sticker that turns red when sterilized.

### Shelf life

5 years



## MATERIALS USED



## PRODUCT RANGE OVERVIEW

| Reference code | Description                     | Size (ml) | Type      | Scaling (ml) | Dead Volume (ml) |
|----------------|---------------------------------|-----------|-----------|--------------|------------------|
| SY3-1SCTU-GIN2 | 3 part syringe 1ml luer slip    | 1         | Centric   | 0,01         | 0,07             |
| SY3-2SC-GIN2   | 3 part syringe 2.5 ml luer lock | 2,5       | Centric   | 0,1          | 0,07             |
| SY3-5SC-GIN2   | 3 part syringe 5ml luer slip    | 5         | Centric   | 0,2          | 0,07             |
| SY3-10SC-GIN2  | 3 part syringe 10ml luer slip   | 10        | Centric   | 0,2          | 0,10             |
| SY3-20ESC-GIN2 | 3 part syringe 20 ml luer slip  | 20        | Eccentric | 1            | 0,15             |
| SY3-20LC-GIN2  | 3 part syringe 20ml luer lock   | 20        | Centric   | 1            | 0,15             |
| SY3-50LC-GIN2  | 3 part syringe 50ml luer lock   | 50        | Centric   | 1            | 0,20             |

## PACKAGING DETAILS

| Reference      | Size | Inner box |        |        |        |             | Outer box |        |        |        |             |
|----------------|------|-----------|--------|--------|--------|-------------|-----------|--------|--------|--------|-------------|
|                | (ml) | Quantity  | L (mm) | W (mm) | H (mm) | Weight (kg) | Quantity  | L (cm) | W (cm) | H (cm) | Weight (kg) |
| SY3-1SCTU-GIN2 | 1    | 100       | 400    | 135    | 65     | 0,40        | 2000      | 550    | 335    | 410    | 8,9         |
| SY3-2SC-GIN2   | 2.5  | 100       | 340    | 120    | 95     | 0,40        | 2000      | 500    | 490    | 360    | 9,8         |
| SY3-5SC-GIN2   | 5    | 100       | 425    | 125    | 115    | 0,68        | 1500      | 590    | 390    | 440    | 12,1        |
| SY3-10SC-GIN2  | 10   | 100       | 450    | 145    | 145    | 1,00        | 1200      | 600    | 445    | 470    | 14,1        |
| SY3-20ESC-GIN2 | 20   | 50        | 415    | 175    | 115    | 0,87        | 500       | 585    | 440    | 425    | 10,2        |
| SY3-20LC-GIN2  | 20   | 50        | 415    | 175    | 115    | 0,87        | 500       | 585    | 360    | 425    | 10,2        |
| SY3-50LC-GIN2  | 50   | 30        | 400    | 183    | 135    | 1,05        | 240       | 550    | 380    | 410    | 10,0        |

### Languages

15 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), Finnish (FI), Polish (PL), Hebrew (HE), and Arabic (AR)

### Transport conditions

Closed and dry

### Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

|                                                                                    | Individual unit | Inner box | Outer box |
|------------------------------------------------------------------------------------|-----------------|-----------|-----------|
| Product name                                                                       | ✓               | ✓         | ✓         |
| Product size                                                                       | ✓               | ✓         | ✓         |
| Quantity (pieces)                                                                  | ✓               | ✓         | ✓         |
| LOT                                                                                | ✓               | ✓         | ✓         |
|   | ✓               | ✓         | ✓         |
| STERILE R                                                                          | ✓               | ✓         | ✓         |
| CE                                                                                 | ✓               | ✓         | ✓         |
|   | ✓               | ✓         | ✓         |
|   | ✓               | ✓         | ✓         |
|   | ✓               | ✓         | ✓         |
|   | ✓               | ✓         | ✓         |
| EC REP                                                                             | ✓               | ✓         | ✓         |
|   | ✓               | ✓         | ✓         |
|  | —               | —         | ✓         |

