



SAFETOUCH™ BLOOD COLLECTION SET

STERILE, SINGLE-USE WITH SAFETY MECHANISM
INTENDED FOR BLOOD COLLECTION



SAFETOUCH™ BLOOD COLLECTION SET

STERILE, SINGLE-USE WITH SAFETY MECHANISM INTENDED FOR BLOOD COLLECTION

- Irreversible, transparent safety mechanism activated by one hand
- Audible click confirms locking
- Does not require a change in puncture technique
- Ultra sharp three-beveled needle to minimize discomfort
- Siliconization of needle minimizes penetration and gliding force
- 20 G luer adapter with self-sheathing needle system
- Soft, large wings for easy grip
- Not made with latex or DEHP
- Available in tubing lengths 19 cm or 30 cm
- Available in 21 G, 22 G, 23 G
- For use by healthcare professionals only

CE 0123

Class IIa Medical Devices

Rule 6

MDD 93/42/EEC

UMDNS 12736

PRODUCT COMPLIANCE

- Complies with the following norms, directives, and regulations:
 - ISO 2859-1
 - ISO 8536-4
 - EN ISO 11135
 - ISO 6009
 - ISO 9626
- Labels contain 12 languages

MANUFACTURING DETAILS

Legal manufacturer: Nipro Corporation

Country of origin: Thailand

STERILIZATION AND SHELF LIFE

Sterilization

EtO (Ethylene oxide)



Single-use only



Indicator sticker


"Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a BLUE sticker that turns RED when sterilized."

Shelf life

5 years



MATERIALS USED



PRODUCT RANGE OVERVIEW

Reference code	Description	Tubing length	Flow rate not less than	Priming volume	Color
BCS-SF-21-19-ET	Blood collection set 21 Gx3/4 (0,80x19 mm)	190 mm	≥ 15,0 ml/min	0,39 ml	Green
BCS-SF-21-30-ET	Blood collection set set 21 Gx3/4 (0,80x19 mm)	300 mm	≥ 15,0 ml/min	0,53 ml	Green
BCS-SF-22-19-ET	Blood collection set 22 Gx3/4 (0,70x19 mm)	190 mm	≥ 12,5 ml/min	0,42 ml	Black
BCS-SF-22-30-ET	Blood collection set 22 Gx3/4 (0,70x19 mm)	300 mm	≥ 12,5 ml/min	0,56 ml	Black
BCS-SF-23-19-ET	Blood collection set 23 Gx3/4 (0,60x19 mm)	190 mm	≥ 7,0 ml/min	0,39 ml	Blue
BCS-SF-23-30-ET	Blood collection set 23 Gx3/4 (0,60x19 mm)	300 mm	≥ 7,0 ml/min	0,54 ml	Blue

PACKAGING DETAILS

Individual unit			Inner box				Outer box					
L (mm)	W (mm)	Weight (g)	Quantity	L (mm)	W (mm)	H (mm)	Weight (g)	Quantity	L (cm)	W (cm)	H (cm)	Weight (kg)
113	110	5,418	50	215	120	102	325,9	500	532	256	244	3,85

Languages

12 languages on all levels of packaging

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), and Finnish (FI)









Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
	✓	✓	✓
	—	—	✓

