



SYRINGE WITH CATHETER TIP

STERILE, SINGLE-USE 50 ML SYRINGE
WITHOUT NEEDLE FOR IRRIGATION



SYRINGE WITH CATHETER TIP

STERILE, SINGLE-USE 50 ML SYRINGE WITHOUT NEEDLE FOR IRRIGATION

- Transparent barrel for precise volume reading
- Smooth plunger motion along entire length of barrel
- Large finger grips for secure handling
- Not made with Latex, BPA, DEHP, or PVC

CE₀₁₂₃

Class Is Medical Devices

Rule 1

MDR 2017/745 Annex IX, chapter I +III

UMDNS 13942

PRODUCT COMPLIANCE

- CE-marked, Class Is Medical Devices, Rule 1, MDR 2017/745 Annex IX, chapter I + III, UMDNS 13942
- Complies with the following norms, directives, and regulations:
 - ISO 7886-1
 - EN ISO 10993-5
 - EN ISO 11607-2
 - ISO 10993-1
 - EN ISO 10993-7
 - EN ISO 13485
 - EN ISO 10993-4
 - EN ISO 11607-1
- Labels contain 15 languages

MANUFACTURING DETAILS

Legal manufacturer: Zhejiang Kindly Medical Devices Co., Ltd.

Country of origin: China

STERILIZATION AND SHELF LIFE

Sterilization

EtO (Ethylene oxide)



Single-use only




Shelf life

5 years



MATERIALS USED



PRODUCT RANGE OVERVIEW

Reference code	Description	Size (ml)	Type
SY3-50CTC-EC	3 part syringe 50 ml catheter tip	50-60	Centric

PACKAGING DETAILS

Reference	Size		Inner box				Outer box				
	L (mm)	W (mm)	Quantity	L (mm)	W (mm)	H (mm)	Quantity	L (cm)	W (cm)	H (cm)	Weight (kg)
SY3-50CTC-EC	260	65	30	260	225	170	360	540	530	485	18.7

Languages

15 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), Finnish (FI), Russian (RU), Turkish (TR), and Polish (PL)









Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
REF	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
	—	—	✓
	✓	✓	✓

