




HUBER NEEDLE INFUSION SET

STERILE, SINGLE-USE, NON-CORING HUBER NEEDLE INFUSION SET



HUBER NEEDLE INFUSION SET

STERILE, SINGLE-USE, NON-CORING HUBER NEEDLE INFUSION SET

- Intended for the infusion of drug solutions into or blood sampling from a reservoir implanted in the body
- Designed for use with a vascular access infusion system
- Non-coring needle with cap
- Needle fixed at 90° angle
- Fixed rotating wings
- Not made with latex or DEHP
- Available with or without Y-site
- Available in 19, 20, 22 G
- Available in 19, 25, 38 mm (¾, 1, 1½ inches, respectively)

CE 0123

Class IIa Medical Devices

Rule 7

MDD 12729

UMDNS 17701

Product Compliance

- Complies with the following norms, directives, and regulations:
 - ISO 10555-3
 - ISO 8536-4
 - ISO 80369-7
 - ISO 7864

Manufacturing DETAILS

Legal manufacturer: Nipro Corporation

Country of origin: Thailand

Sterilization and Shelf Life

Sterilization

EtO (Ethylene oxide)



Single-use only



Indicator sticker

Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a violet sticker that turns to blue when sterilized.


Shelf life

3 years



Materials Used

**Not made with
Latex or DEHP**




Lubricant

Material: Silicone

Glue

Material: Epoxy



PRODUCT RANGE OVERVIEW

Reference code	Description	Needle outer diameter	Needle length	Tube length volume	Y-site
HCNS-1919-ET	Huber coreless Ndl Set 19 G 19 mm	19 G	19 mm	100 mm	No
HCNS-1925-ET	Huber coreless Ndl Set 19 G 25 mm	19 G	25 mm	100 mm	No
HCNS-2019-ET	Huber coreless Ndl Set 20 G 19 mm	20 G	19 mm	100 mm	No
HCNS-2025-ET	Huber coreless Ndl Set 20 G 25 mm	20 G	25 mm	100 mm	No
HCNS-2038-ET	Huber coreless Ndl Set 20 G 38 mm	20 G	38 mm	100 mm	No
HCNS-2219-ET	Huber coreless Ndl Set 22 G 19 mm	22 G	19 mm	100 mm	No
HCNS-2225-ET	Huber coreless Ndl Set 22 G 25 mm	22 G	25 mm	100 mm	No
HCNS-1919-Y-ET	Huber coreless Ndl Set 19 G 19 mm Y site	19 G	19 mm	200 mm	Yes
HCNS-1925-Y-ET	Huber coreless Ndl Set 19 G 25 mm Y site	19 G	25 mm	200 mm	Yes
HCNS-2019-Y-ET	Huber coreless Ndl Set 20 G 19 mm Y site	20 G	19 mm	200 mm	Yes
HCNS-2025-Y-ET	Huber coreless Ndl Set 20 G 25 mm Y site	20 G	25 mm	200 mm	Yes
HCNS-2038-Y-ET	Huber coreless Ndl Set 20 G 38 mm Y site	20 G	38 mm	200 mm	Yes
HCNS-2219-Y-ET	Huber coreless Ndl Set 22 G 19 mm Y site	22 G	19 mm	200 mm	Yes
HCNS-2225-Y-ET	Huber coreless Ndl Set 22 G 25 mm Y site	22 G	25 mm	200 mm	Yes
HCNS-2238-Y-ET	Huber coreless Ndl Set 22 G 38 mm Y site	22 G	38 mm	200 mm	Yes

Packaging Details

Without Y-site

Individual unit		Inner box				Intermediate box				Outer box				
L (mm)	H (mm)	Quantity (pcs)	L (mm)	W (mm)	H (mm)	Quantity (pcs)	L (mm)	W (mm)	H (mm)	Quantity (pcs)	L (mm)	W (mm)	H (mm)	Weight (kg)
160	100	25	160	115	102	50	235	162	105	1000	558	498	361	10,8

With Y-site

Individual unit		Inner box				Intermediate box				Outer box				
L (mm)	H (mm)	Quantity (pcs)	L (mm)	W (mm)	H (mm)	Quantity (pcs)	L (mm)	W (mm)	H (mm)	Quantity (pcs)	L (mm)	W (mm)	H (mm)	Weight (kg)
160	100	25	160	115	102	50	235	162	105	1000	558	498	361	10,8

Languages

26 languages of the IFU available:

(BG) Bulgarian, (CS) Czech, (DA) Danish, (DE) German, (EL) Greek, (EN) English, (ES) Spanish, (ET) Estonian, (FI) Finnish, (FR) French, (HR) Croatian, (HU) Hungarian, (IT) Italian, (LT) Lithuanian, (LV) Latvian, (NL) Dutch, (NO) Norwegian, (PL) Polish, (PT) Portuguese, (RO) Romanian, (RU) Russian, (SK) Slovak, (SL) Slovenian, (SR) Serbian, (SV) Swedish, (TR) Turkish.

Transport and storage conditions

Keep dry and away from sunlight.

Label details

	Individual unit	Inner box	Intermediate box	Outer box
Product name	✓	✓	✓	✓
Product size	✓	✓	✓	✓
Quantity (pcs)	—	✓	✓	✓
REF	✓	✓	✓	✓
LOT	✓	✓	✓	✓
	✓	✓	✓	✓
STERILE EO	✓	✓	✓	✓
CE	✓	✓	✓	✓
	✓	✓	✓	✓
	✓	✓	✓	✓
	✓	✓	✓	✓

	Individual unit	Inner box	Intermediate box	Outer box
	✓	✓	✓	✓
EC REP	✓	✓	✓	✓
	—	—	—	✓
	—	—	—	✓
	—	—	—	✓
	—	—	—	—
	—	—	—	✓
	✓	✓	✓	✓
Not made with LATEX	✓	✓	✓	✓
Not made with DEHP	✓	✓	✓	✓

