



LUER ADAPTER

STERILE, SINGLE-USE MALE LUER ADAPTER
FOR BLOOD COLLECTION



LUER ADAPTER

STERILE, SINGLE-USE MALE LUER ADAPTER FOR BLOOD COLLECTION

- Multiple draws with a single vein puncture to optimize efficiency
- Standardized luer connection for universal fit
- Not made with latex, DEHP or PVC
- Available in 20 G
- For use by healthcare professionals only

CE 0123

Class IIa Medical Devices
Rule 2
MDD 93/42/EEC
UMDNS 12736

PRODUCT COMPLIANCE

- Complies with the following norms, directives, and regulations:
 - BS EN 556-1
 - BS EN 1041
 - ISO 7864
 - ISO 9626
 - ISO 10993-1
 - ISO 10993-7
 - ISO 10993-10
 - ISO 10993-11
 - ISO 11135
 - ISO 11607-1
 - ISO 11607-2
 - EN ISO 11737-1
 - EN ISO 13485
 - ISO 15223-1
 - ISO 14644-1
 - ISO 14644-2
 - ISO 14698-1
 - EN ISO 14971
 - EN 62366
 - ISO 80369-7

MANUFACTURING DETAILS

Legal manufacturer: Nipro Corporation

Country of origin: JAPAN

STERILIZATION AND SHELF LIFE

Sterilization

ETO (Ethylene oxide)



Single-use only



Indicator sticker

Each outer box (or shipping carton) contains a chemical indicator that ensures sterility. This indicator comes in the form of a [pink] sticker that turns [blue] when sterilized.

Shelf life

3 years



MATERIALS USED



Not made with latex, DEHP or PVC

PRODUCT RANGE OVERVIEW

Reference code	Description	Container	Cap	Label
LM-20-EJ	Luer adapter 20G	✓	✓	✓

PACKAGING DETAILS

Individual				INNER BOX					OUTER BOX (mm)				
L (cm)	W (cm)	H (cm)	Weight (g)	Quantity	L (cm)	W (cm)	H (cm)	Weight (kg)	Quantity (pcs)	L (cm)	W (cm)	H (cm)	Weight (kg)
0,7	7,5	0,7	1,96	100	8,1	16,4	5,2	0,2	10.000	90	44,5	24,5	24,6

Languages

12 languages

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Finnish (FI), Swedish (SV), Danish (DA), Norwegian (NO), Greek (EL)


Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Intermediate box	Outer box
Product name	✓	✓	✓	✓
Product size	✓	✓	✓	✓
Quantity (pieces)	—	✓	✓	✓
REF	✓	✓	✓	✓
LOT	✓	✓	✓	✓
	✓	✓	✓	✓
STERILE EO	—	✓	✓	✓
CE	—	✓	✓	✓
	—	✓	✓	✓
	—	✓	✓	✓
	—	✓	✓	✓
EC REP	—	✓	✓	✓
	—	✓	✓	✓
	—	—	✓	✓

