



SYRINGE WITH LUER LOCK, WITHOUT NEEDLE

STERILE, SINGLE-USE SYRINGE FOR INJECTION OR ASPIRATION

Manufactured in Guangdong, China



SYRINGE WITH LUER LOCK, WITHOUT NEEDLE

STERILE, SINGLE-USE SYRINGE FOR INJECTION OR ASPIRATION

- Transparent barrel for easy recognition of air bubbles
- Smooth plunger motion along entire length of barrel
- Large finger grips for secure handling
- Latex-free, BPA-free, DEHP-free, PVC-free
- Available in centric tip

CE₀₁₂₃

Class I Medical Devices

Rule 1

MDD 93/42/EEC

UMDNS 13929

PRODUCT COMPLIANCE

- CE-marked, Class I Medical Devices, Rule 1, MDD 93/42/EEC, UMDNS 13929
- Complies with the following norms, directives, and regulations:
 - ISO 594-1
 - EN 1707
 - ISO 7886-1
 - ISO 13485
 - ISO 10993-1
 - ISO 10993-4
 - ISO 10993-5
 - ISO 10993-7
 - ISO 11607-1
 - ISO 11607-2
 - EN ISO 13485
 - ISO 15223-1
- Labels contain 15 languages

MANUFACTURING DETAILS

Legal manufacturer: Guangdong Kindly Medical Device Group Co., Ltd.

Country of origin: China

STERILIZATION AND SHELF LIFE

Sterilization

EtO (Ethylene oxide)



Single-use only




Shelf life

5 years



MATERIALS USED



DEHP-free
BPA-free
Latex-free
PVC-free

PRODUCT RANGE OVERVIEW

Reference code	Description	Size (ml)	Tip type
SY3-3LC-GEC	3 part syringe 3 ml luer lock G	3	CENTRIC
SY3-5LC-GEC	3 part syringe 5 ml luer lock G	5	CENTRIC
SY3-10LC-GEC	3 part syringe 10 ml luer lock G	10	CENTRIC
SY3-20LC-GEC	3 part syringe 20 ml luer lock G	20	CENTRIC
SY3-30LC-GEC	3 part syringe 30 ml luer lock G	30	CENTRIC
SY3-50LC-GEC	3 part syringe 50 ml luer lock G	50	CENTRIC

PACKAGING DETAILS

Reference	Size	Inner box				Outer box				Net weight	Gross weight
	(ml)	Quantity	L (mm)	W (mm)	H (mm)	Quantity	L (cm)	W (cm)	H (cm)	(kg)	(kg)
SY3-3LC-GEC	3	100	420	125	70	3200	58	44	53	13,3	18,7
SY3-5LC-GEC	5	100	420	125	93	2400	58	44	46	14,0	18,8
SY3-10LC-GEC	10	100	420	165	110	1500	58	44	46	13,7	17,8
SY3-20LC-GEC	20	50	215	210	140	800	58	44	46	11,9	15,6
SY3-30LC-GEC	30	50	285	210	153	600	58	44	50	12,0	15,8
SY3-50LC-GEC	50	30	280	210	140	360	58	44	46	12,2	15,6

Languages

15 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), Finnish (FI), Russian (RU), Turkish (TR), and Polish (PL)





Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
REF	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
	—	—	✓
	✓	✓	✓

